

Product datasheet for **MC203563**

Saa3 (NM_011315) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Saa3 (NM_011315) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Saa3
Synonyms:	AV098916; I7R3; Saa-3
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC055885 TGGGGTCCCAGCAGGAGCTCGCAGCAGCAGCAGGATGAAGCCTTCCATTGCCATCATTCTTTGCATCTTG ATCCTGGGAGTTGACAGCCAAAGATGGGTCCAGTTCATGAAAGAAGCTGGTCAAGGGTCTAGAGACATGT GGCGAGCCTACTCTGACATGAAGAAAGCTAACTGGAAAACTCAGACAAATACTTCCATGCTCGGGGGAA CTATGATGCTGCCCGGAGGGGTCCCGGGGAGCCTGGGCTGCTAAAGTCATCAGTGATGCCAGAGAGGCT GTTCAGAAGTTCACGGGACATGGAGCAGAGGACTCAAGAGCTGACCAGTTTGCCAACGAGTGGGGCCGGA GTGGCAAAGACCCCAACCACTTCCGACCTGCTGGCCTGCCTAAAAGATACTGAGTTTTCTTCTCCTGTTG TTCCCTGTGATGCTGCCCCCGAGAAGAGGGGCAACTACTGGGTTGAGATATTTTCTAAATCTGGATCC CTAAACATCCCAATGTGCTGAATAAATACTTGTGAAATGCAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_011315
Insert Size:	369 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC055885](#), [AAH55885](#)

RefSeq Size: 555 bp

RefSeq ORF: 369 bp

Locus ID: 20210

UniProt ID: [P04918](#)

Cytogenetics: 7 30.51 cM

Gene Summary: Major acute phase reactant. Apolipoprotein of the HDL complex.[UniProtKB/Swiss-Prot Function]