

## Product datasheet for **MC203624**

### Nifk (NM\_026472) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Nifk (NM_026472) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Nifk
Synonyms:	A1852665; C130020J04Rik; Mki67ip
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC058559

```

CGTGGATTGACTTGTTCGGGTTCTCCAGCATGGCTGGGTTAGCAGGCCCGCGAAGCCGAGTCTGGCATT
GAACCCGCAGGAAGACTCGCAGTTTGAAGAGGCTCTGACGCAGATCCAAGGGCGCACGAAGAAGCCCTCAG
CAAAAGAAGAAAGAAAACTTAACCGTGGAGTAGTTTATCTGGGCCACCTACCTTCAACACTTTCTGAAA
GCCACATCTATAACTACTGTGCCAGTTTGGCGCATTAGCAGATTGAGATTGTCTAGAAGTAAGCGGAC
TGGAAACAGCAAAGGCTATGCCTTTGTGGAGTTTGTGAGTCTGAGGATGTTGCCAAGATAGTTGCAGAAACA
ATGGATAACTACCTTTTTGGTGAAGACTTCTATCATGTAATTTATGCCACGAAAAAAGTCCATAAAG
ACCTTTTCAGCCAGAGGAATGCTCTGTTTCACCGGCCATCATTCCAGCAGTGAAGCGCTATAATCGGAA
ACGAGGACATTTGCAGATGTTGAAGATGGAGTATAGGTTTAAGAAGAAGGAAAAATTAATCCGGAAGAAA
CTAGCTCAAAGGGGATTGATTATAGTTTTCCCTCACTGGTCTACCTAAGCCAAAGAACATTGCTGTCCG
CTCATAGAGATTCTGAGGTAACCAGTTTTACCTGACCAAAAAGAAGGGCTTTCCGGGTGAACTAGGAG
AAAGGAGAAGATGATGAAAGAGGACATTTCCAATAACATTCTCAAGAAGCGTAAGAGGTCTAGGAGGAAG
AAGTCGTCTGTGGACAGCCAGGGTCCCACACCAGTTTGTACACCGACGTTTTTGGAGAGGCGAAAAATCGC
AGGTGATGGAAGTCGGTGGCGATAAAGATGATGAAATCATTCTCAAGCTGCCTGTGCCCCAGTGAAGA
AGACACACAAAAGACCCCAACATCTGCAAGCCCTGGGGGAAAAAGACCGAGAAAAAGGAAGCAAGCAG
TGATCCCAGCAGCGTTCTGTGCATAGTTCTAGGCAGATGTGTTGTTATAAGGGACTACTATTTTACTA
AGTGCAGTGACTTGCAAGCACAGACAAGACTGTTGGAGGAAGAATACATGTTTGGTTTTCCAGTGAGGAA
AGCCACTGGGAACCTGACTGTGCCCTTTCTGTCTGAATTTGTCTTCAGGCTTCATCGGTCCCTCAGTT
TTACTATCCTTATGAAACTCAGGCTTAACCTACCACCAACAACCGATTGGTGTGTTTGGACGGTGTGTA
GTTGTGTATCTAAACCCCTCCAAAATATGTCAGGCTTCAGCCTCAGACATCTGAAACTGGAATGAAATC
CAGAGAGTGAGAGGGCAGTTCTGGTTTAGATTTCCACATGTGTCCAGGAGCTGAGGAGGAGTCCGAATG
ACCCACAGTTACTGCAGCTAAGCTAGTGTGTAGGTGTCACAAATCAGTTTCAGCAGTTGTATTCAAGATT
CACTGACCAGACTTTTCTACAGTTTATTTGTGTACAAGAGAACCTAAGAATTTGTGGGGAAGTAAGAAA
TGTCAACATCTCCTGTTGAAGTTTCTGTACCCTTTCTATGCTTTTGTTTTACTAATTTGGGTTTCTTGA
ATATAACTGTCAATCTCCTTGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAA
AAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_026472
Insert Size:	954 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
RefSeq:	<a href="#">BC058559</a> , <a href="#">AAH58559</a>
RefSeq Size:	1720 bp
RefSeq ORF:	954 bp
Locus ID:	67949
UniProt ID:	<a href="#">Q91VE6</a>
Cytogenetics:	1 E2.3