

Product datasheet for **MC203737**

Katna1 (NM_011835) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Katna1 (NM_011835) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Katna1
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC009136

```
CGACCCGGAGCTGCGCGGGTAGCCTGGGCTGCGCTAGCTTACATATGACCATGAGTCTTCAAATGATTG
TTGAGAATGTAAATTTGGCTCGTGAATATGCACTGCTGGGAACTATGACTCTGCAATGGTCTACTATCA
GGGAGTTCTTGACCAAATGAACAAGTACCTGTACTCAGTCAAAGATACACACCTCCGTCAGAAATGGCAA
CAGGTTTGGCAGGAAATAAATGTGGAAGCTAAGCAAGTTAAGGATATCATGAAAACTCGAGAGCTTCA
AACTAGACATCACTTCTTTGCAAGCTGCACAACATGAACCTCCAGCTGCTGAAGGAGAAGTCTGGTCCTT
GCCTGTACCTGTTGAGCGAAGACCCTTACCAGGACCTAGGAAACGCCAGTCTTCTCAGCACAGTGACCTT
AAGCCACACAGTAACCGGCCAAGCACAGTCGTGAGAGCTCATCGGCCGTCTCCACAGAATCTTCACAATG
ACAGAGGGAAAGCTGTTTCGTAGTCGTGAAAAGAAAAGAACAGAGTAAGGGAAGAGAGGAAAAGAACAAAT
ACCTGCTGCAGTTACAGAACAGAGGCAAATAAGTTTGATGGCACTGGGTATGACAAGGACTTAGTAGAA
GCTTTGGAAAGAGATATCATTTCTCAGAATCCCAACGTTTCGGTGGTATGATATTGCTGATTTAGTAGAAG
CCAAAAAGTTGCTTCAGGAAGCAGTGGTGTACCAATGTGGATGCCAGAATTCCTTAAGGGCATTAGGAG
ACCATGGAAGGTGTGCTGATGGTTGGCCACCTGGCACTGGAAGACCTTCTAGCTAAAGCAGTTGCC
ACGGAATGCAAGACAACGTTCTTCAATGTCTTTCATCAACTCTGACGTCCAAATACAGAGGAGAGTCTG
AGAAGCTTGTCCGTCTTCTGTTTGAATGGCTCGGTTTTACTCCCAGCTACCATATTTATTGATGAGAT
AGATTCCATTTGCAGCCGCCGAGGAACCTCTGAAGAGCATGAAGCCAGCCGAAGGATGAAAGCTGAGCTG
CTGGTTCAAATGGATGGTGTGGAGGTGCCTCAGAGAACGATGACCTTCCAAATGGTTATGGTTCTGG
CGGCTACTAATTTCCCATGGGACATAGATGAGGCTTTAAGGCGTCGTCTGGAGAAAAGAATCTATATTCC
ATTACCATCAGCAAAAGGCAGGGAAGAACTGTTGAGAATAAGTCTGCGAGAGCTGGAATTTGGCTGATGAC
GTTAACCTTGCAGCATAGCAGAGAACATGGAAGGCTACTCGGGTGCAGACATTACCAATGTATGCAGGG
ACGCATCGTTGATGGCCATGAGACGGCGTATTGAGGGTCTAACCCAGAAAGAAATCCGGAACCTTTCAAG
AGAAGCAATGCACATGCCTACCACCATGGAGGATTTTGAATGGCTTTAAAGAAGATTTCTAAATCAGTA
TCAGCTGCAGATATTGAAAGATACGAGAAATGGATAGTTGAGTTTGGATCATGCTAGCTATGTTCTCACA
TGTAAGCTGTGAGGAATCTGCCTTAAGTCTTTGAAAAATAATTAATGTAGAATTCAGTGTACGGAGTG
CTATATTTTTTAACTTTCATAATGGTGTGATTTGGAAATCTTTCTGGTTCTGAATAAAGACTATTTT
TAAGCTGCAAAAAAAAAAAAAA
```

Restriction Sites: RsrII-NotI


[View online »](#)

ACCN:	NM_011835
Insert Size:	1476 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC009136, AAH09136</u>
RefSeq Size:	1703 bp
RefSeq ORF:	1476 bp
Locus ID:	23924
Cytogenetics:	10 A1