

Product datasheet for **MC203743**

Lenep (NM_020517) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lenep (NM_020517) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Lenep
Synonyms:	AU018355; Lep503
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC006806

```

CTGAAAGGCTGAGTCAGCCTGCTCAGAGGCCAGAACCTTCTTTTGTAGACCCTAGAGCCTTTGTTTTT
AGTCATGAGGCCCTGAAGGGGTGGGACAGTGAATAAAAGGTCCCAAAGGCAGTGCCACACCCTTCTGCT
GAGCTGCTGTCCACATGAGGCCACCAACGCAGCCCTGACCCAGGCCCTGCCCTTCTCTCAGAGATGC
CCTTCGAGGCACTGGCCTCCAGGTGCCTGTCAATAAAATGGGCACAGGTTGGGAGGGCATGTATCGAAC
CTGAAGGAAGTCGCCTACATCCTTCTTTGTTGCTGGTGTATCAAGGAACTGCTGGATTAATGGTGGCAGG
GAACACCTCCTGTCAACCCTGACACCAAGACCTGCATCCCTCAGGTGGGAGGGGAGGGCAGACTCTCAG
GAGGTGTGGTGCAGCTTCTGGGTGGAACCTGCAGCCAGCTCTCCCATCTTTTCTATAGCCCTAAGAC
TCCTGGAGCCAGCTGGCCTTACCAAGGATGGGTCAAAAGAACTGAAGAAAGAACCACGAACCTTGCTC
TGAGCTGTCCATCCATCTGACCTGCCTCACCCTGGGAATCTCTCCTCAGCTCCCATCTCTGCTGTGCTC
CTGGTACGCATGAGCTTGGCTCAGACTCCATTGAGACTACGGCAGGATCTAATATTCCTCAAAAATTGG
CACAGTGAGCCACCAACCCACGTACTTGAAGAACTAGGAATGGACTGCCCTGAAATGCCTGGTAAATCA
CTCCAATAATAATAATAATAAGGCCTCGAGGCCAAAAATTCTGAACATTCAAAGGGGAATGAGGGCAT
GTGTGCTTGAAAGGGCTTAGAGCCAGAGAGCAGGAAGAAGGACAGATCTCGGGAGTAGTGAGACCAAGA
ACCAGCTTTGCTGAATCTCTGTTGCCAGGAGGGGAAGAGCAAAACAGAAAAGCCATGTTTCTCCCACTTT
TGTGCAGAGAGGACAGAAATGGAGCAAGCCCTATTCTATGAGGCCCTGGGATGCTAGGATGAGGACCCAAA
ACTCAGCCTTTTCCAGAAGATGTGACCTAGGCAGCTCCCTGTCTAAGGAACACGGGAAGTGTACTCAG
CTAGTCCCAGTTCAGATACTAAGACACCCTAAGAGAAAACACCCATGAGGTCATACCACCCCTTCAAG
CCTGTACAGCACACGTGTGACAGCACATAATCTGCGTGAGTGAATGGCCTCCTGACTTGGCGTTGGCAG
CCACCAGCACTTGTCTTTCTTTGAACTGTTAGCAGCCTGTGCACATGCGTATACATCAACTGTTTCATGCG
CCATCCCCACCCACCCCATCCCCACCCTTACCCTTAGTAGCACTCTTTGGCCTAGCCATTGCCACCT
CTTACCTGCAGCACTGACCCGTAGGTAGTATATAGCTCTCTGCCCATCTCCTGGAGCCAGCTCTTTTTT
GGCCTTTTGTAGTTAGTCTCTTTTTCTCGGCCTGTTTTATTGTTGGTTATCTGGTTGGTTGGTTG
GTTGGTTGAGTTTCCGTGTTGGGACTGAACCCAGGACTTTGTACCTACTCCTCTCCCTATACCACAGCC
TGATTTTTGTTGCATACAGATACTATATGTAAGTTGTTGGTAACCCACACCTTAATTCAGCACTTG
GGAGGCGGAGGCTGGAAGATTTGAGACCATCCTGATCTATACAGCACATTCTTCAGTGAATAAGTGAAT
ATGAATTGGATAATCAAATGTTAAAATTAGTGAAGAAGCAAAAAAAAAAAAAAAAAAGGGCGGCCG

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_020517
Insert Size:	186 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC006806</u> , <u>AAH06806</u>
RefSeq Size:	1815 bp
RefSeq ORF:	186 bp
Locus ID:	57275
UniProt ID:	<u>Q9WVB6</u>
Cytogenetics:	3 39.1 cM