

Product datasheet for MC203758

Msmo1 (NM_025436) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Msmo1 (NM_025436) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Msmo1
Synonyms:	150001G16Rik; C78600; DESP4; ERG25; Sc4mol
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC006802

```
GCCAGGACTTGCCTTAGACCGCTGTCTGTGCGGGGACGCCAGGAATAAAAGGTCCTCGCAGAGATT
GAAAATGGCCACAAACAAAAGTGTGGCGTGTTCAGCTCTGCATCCCTGGCTGTGGAGTATGTGGACTCA
CTGTTGCCGGAGAATCCTCTGCAGGAGCCCTTTAAGAATGCTTGGGTCTACATGTTGGACAACATACAA
AGTCCAGATCGCAACCTGGGGCTCCCTGATAGTTCACGAAGCCATCTATTTCTTGTCTCTTTACCTGG
ATTCTATTCCAGTTTATACCTTACATGAGAAAGTACAAGATTCAAAGGATAAACCAGAAAACCTTCGAA
GGCCAGTGGAAAGTGTGAAAAAATTCTATTTAATCATTTTTTTATCCAGTTGCCTCTGATTTGTGGAA
CATACTATTTTACAGAGTTTTTCAATATTCCTTATGATTGGGAAAGAATGCCAAGATGGTACTTAACTTT
GGCAAGGTGTTTGGGCTGTGCAGTCATTGAGGACACCTGGCATTATTTCTGCACAGACTCCTTCACCAC
AAGAGAATTTATAAATACATTCAAAAGTCCATCAGAGTTTCAGGCTCCATTTGGAATCGAAGCAGAAT
ACGCACATCCCTTGAAACCTGATTCTTGGGACAGGATTTTTTCATCGGAATTGTCTTTGTGTGATCA
CGTAATTCCTGTGGGCATGGGTGACCATACGTTTGTGGAACCATCGATGTGCACAGCGGTTATGAT
ATTCCTCTCAACCCGCTGAACTTGGTTCCTTTTTATACGGGCGCTCGGCACCATGATTTCCACCACATGA
ACTTCATTGGGAATTATGCCTCCACCTTACGTGGTGGGATAAATCTTTGGGACAGATGCCAGATATCA
TGCTTACATTGAAAAATCGAAGAAGCTTGGGAAAAAGTCTGATTGAGCACCCCACTCGACGCCTTCCTC
AACTGCAGCTGGCCACTAACCTGCTTCTGGGACTTGGCTGCCAAGTGATAAAACGTTAACGACCTTT
AAGTACTCTTCCAGATGAGAGCTCTTCTACCTTACGTCCAAGATCTGTATATGTGGGAACAAATAAGTAA
ATATATACATATATATATATTTAAGTACAGTTTTTATCAGGAGTTTTAAAGGCCATATTGCTAACCT
CACAAAAGGTTTCGGGAAGTGGGAAAGTATTAATGTACAACCTCACCCTTCAAGGAGGCTGA
AGTTGACGTTGCCTTTTAAAGCCTTTTACATACACTGGTCAAGTTCGAAAAATTCTCAACAATAATGTTTGC
CTTCTAGTTTAAAGACATGGTGTATTTTAAAGAAATACCTGAAGTGAAGCTATCACGAAGTGAAGTGTATCAC
CTGGTGAAGCAGTACTGAAACCTGCTCTGGTCTGGGTTTACCAGGTTAGGACTATTCTTTCTGTGAAG
GATTGTTTATTGATAGTAACTCTATTGGCATCTTGTAATTTGGCTCATTTCCTACGTGATGTACATAGAA
TATTATAGCCTAAAGATTTTCCAGTTACTTGGCGTCTTTTTAAATTTGAGGAGCCCACTTCCACTGTCCA
CCCCTTGTCTTTGTGTGCAGGTGATCATTGAAATCGATCTTGGAGCTTGAATTCACCTTCTGTCTGTT
TCCGGAGCTGGTAAATAAATGTGATTTTTGTCCGATCCTTCTGGTTTTTATGAAGATTAGCTGCCATGT
TTTATTCTTATAAGGAAAAGCACAGAGATCATTAAAGTAAGTTTTACAATAAAAAAAAAAAAAAAAAAAAAA AAAAAA
```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_025436
Insert Size:	882 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC006802 , AAH06802
RefSeq Size:	1827 bp
RefSeq ORF:	882 bp
Locus ID:	66234
UniProt ID:	Q9CRA4
Cytogenetics:	8 B3.1
Gene Summary:	Catalyzes the first step in the removal of the two C-4 methyl groups of 4,4-dimethylzymosterol.[UniProtKB/Swiss-Prot Function]