

## Product datasheet for MC203766

### Gnb1 (NM\_008142) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Gnb1 (NM_008142) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gnb1
Synonyms:	AA409223; C77571; Gnb-1
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC013058

```

GCCGCGCCCATCCTGATGGCTGGAGGCGTCTGAGGGGCGGACGGCGACGGCGCGCGGGCGGGCGGGCGGTG
GCGGGCGCGGGCGGCGAGGGCAGCGGGGCTGGGACTGGAACAGCCGCGCGGACCCGACCAGCCTCGC
CGACTCCTGCCCCAGGCCGCGAGCGCGACGACCAGGCGCGCGGCTGGGCGCCAGACCGGGCGAGGGCGA
GATAGACTTACATGTGTGGGAAACCAGACCAGAAACACTTCTGAGTCAAGAGCTTGGATTCTGAAGTGTG
GAAAGCACTGAGACGTGAAGATGAGTGAACCTTGACCAGCTGCGGCAGGAGGCCGAGCAACTGAAGAACCA
AATTAGAGATGCTCGTAAAGCGTGTGCCGATGCGACTCTTTCTCAGATCACAACAATATTGATCCAGTG
GGAAGAATCCAAATGCGGACCAGGAGAACAAGTGGGGGCGATCTGGCAAAGATTTATGCCATGCACTGGG
GCACAGACTCAAGGCTCCTTGTGACGCGCTCTCAGGATGGAAAACATCATCTGGGACAGTTATACCAC
AAACAAGGTTTCATGCCATCCCTCTGCGCTCCTCTGGGTCATGACCTGCGCATAACGCTCCTTCTGGGAAT
TATGTGGCCTGTGGTGGCCTGGATAACATCTGCTCCATTTACAACCTGAAAACCTGTAAGGGAATGTGC
GTGTGAGTCGTGAGCTGGCGGGACACAGGTTATCTGTCTGTTGCCGTTTCTGGATGACAATCAGAT
AGTTACCAGTTCTGGAGACACCACATGTGCCCTGTGGGACATCGAGACTGGCCAGCAGACAACCACATTT
ACTGGACACACTGGAGATGTCATGAGCCTGTCTCTTGCTCCTGACACCAGACTGTTTGTCTCTGGTGCTT
GTGATGCTTCAGCCAAGCTCTGGGATGTCCGAGAAGGGATGTGCCGGCAGACCTTTACAGGACACGAGTC
TGACATCAATGCCATATGTTTCTTTCCCAATGGCAATGCCTTTGCCACTGGCTCAGACGATGCCACATGC
AGGCTGTTTGACCTCCGTGCAGACCAGGAGCTCATGACCTACTCCCATGACAACATTATCTGTGGTATCA
CATCTGTTTCTTCTCCAAGAGTGGCCGCTCCTCCTTGTGGTATGATGACTTCAACTGTAATGTCTG
GGATGCACTCAAAGCTGACAGAGCAGGTGCTTAGCTGGACACGACAACCGAGTCAGCTGCTTGGGGGTG
ACTGATGATGGCATGGCTGTGGCAACAGGGTCCTGGGACAGCTTCTCAAGATCTGGAACATAAAGCCAGG
AGCAGGTGGATGGATGCCATGGCGACTGGAAGACCATTCCAACCTTGAATCGTTACAGCGATAGCATATC
CTATCCAACCATACTAACGTGGACACCCACACCGCCCTCAGAACTTCAAAGGGCAAGATTTTTTTTCC
TTCACTTATTGCTGAAACCAAGAGCACAATTCACATTAAGAGAAGGATTTTTGTGCTGTAATAACTAAACA
ATTGTGCATTTTTTTGGGGCATTGTTTTTGGTTTTATTTTTGTTATTTTTGTTTTGTTTTGTTTTGTT
TTTTTGGTTTCATTTTTGTTTTGGTTGGTCCAATTTTAAAGGAAATACTACAATAAAAATGTTAACCA
AAGGTA
  
```



[View online »](#)

<b>Restriction Sites:</b>	RsrII-NotI
<b>ACCN:</b>	NM_008142
<b>Insert Size:</b>	1023 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC013058</a> , <a href="#">AAH13058</a>
<b>RefSeq Size:</b>	1714 bp
<b>RefSeq ORF:</b>	1023 bp
<b>Locus ID:</b>	14688
<b>UniProt ID:</b>	<a href="#">P62874</a>
<b>Cytogenetics:</b>	4 86.17 cM
<b>Gene Summary:</b>	<p>Guanine nucleotide-binding proteins (G proteins) are involved as a modulator or transducer in various transmembrane signaling systems. The beta and gamma chains are required for the GTPase activity, for replacement of GDP by GTP, and for G protein-effector interaction (By similarity).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript. Variants 1-3 encode the same the protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>