

Product datasheet for **MC203789**

Ola1 (NM_025942) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ola1 (NM_025942) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ola1
Synonyms:	2510025G09Rik; 2810405J23Rik; 2810409H07Rik; Gtpbp9; PTD004
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC011318

```

GGCTGCCGCGACACCGGTCTTTGCAGTCCGGGACCCCGTGCGATCGTTAGCCCGGTGCACGATGCCCCC
GAAAAAGGGAGGAGATGGAATTAACCGCTCCAATTATTGGAAGATTTGGAACCTCACTGAAAAATTGGT
ATCGTTGGATTGCCAAATGTTGGGAAGTCTACCTTCTTCAATGTATTAACCAATAGTCAGGCTTCAGCAG
AAAACCTCCCATTTCTGCACTATTGATCCTAATGAGAGCAGAGTGCCTGTGCCAGATGAGAGGTTCTGACTT
TCTTTGCCAGTACCATAAACAGCAAGCAAGATTCTGCTTTCTCTAAATGTAGTGATATTGCTGGCCTT
GTAAAAGGAGCTCACAATGGGCAGGCTTGGGGAATGCCTTTTTATCTCATATTAGTGCCTGTGATGGGA
TTTTTCATCTAACACGTGCTTTTGAAGACGATGATACACGCATGTTGAAGGAAGTAGATCCCATCCG
AGACATAGAGATAATTCACGAGGAGCTCCAGCTGAAGGACGAGGAGATGATTGGGCCCATCTAGACAAA
CTAGAGAAGGTGGCCGTGAGGGGAGGCGACAAGAACTGAAGCCGAGTATGACATAATGTGCAAAGTCA
AGTCTTGGGTTATAGATCAAAAGAAACCTGTTCTGCTTCTATCATGACTGGAATGACAAAGAGATTGAAGT
GTTGAATAAACACTTATTTTTGACTTCAAACCAATGGTCTACTTGGTGAATCTTTCTGAAAAAGACTAT
ATAAGAAAGAAAAACAAATGGTTGATAAAAAATTAAGAGTGGGTGGACAAGTATGACCCAGGCGCCTTGG
TCATTCCTTTTAGTGGGGCCTTGGAACTCAAGTTGCAAGAATTGAGTGTGAGGAGAGACAGAAGTATCT
GGAAGCGAACATGACACAAAGTGCTTTGCCAAAGATCATTAAAGCTGGATTTGCAGCACTCCAAC TAGAA
TACTTTTTCACTGCAGGCCCAGATGAAGTGCCTGCGTGGACCATCAGGAAAGGAACGAAGGCTCCTCAGG
CTGCAGGAAAGATTACACAGATTTGAAAAGGGATTCAATATGGCTGAAGTAATGAAGTATGAAGACTT
TAAAGACGAGGGTTCTGAAAACGCAGTCAAGGCTGCTGGAAAAACAGACAACAAGGCAGGAATTATATT
GTTGAAGATGGAGATATTATCTTCTTCAAATTTAATACACCTCAGCAACCAAAAAAGAAATAAATTCAGT
TATTGCTAAGATCGACGTACAACCTCTAAAATGCATATATTGGGGCTTGGGGGGATGGGTGTTTTTAAA
TTAAAGTTTCTGAAATCTAAAGAGACAAAGTTAGGAATGGGAATCTTCTACAAACAAATTTTTTTATTT
GTTTTAAATTTGAAAATACTATGTACCTTCAATGAAATGCAAGTTACTAAGCGTGAAGAGCTTTGCTTT
CTCACGTGATTAAGACCCTACTGCACACTGTAGAGGCTTTTCAAGCCATATTACTCTCCTGATGCTTCAT
TAATCTCCATCATGTATGCCAGGCTGACACGTTTGACAGTGAAGACAATGTGGCTTGCTCCTTTTTGAA
TCTAAAGATAATGCATGTTTTACAGTACTCCAGATGTCTACACTCAATAAAACATTTGACAAAACCAAAA AAAAAAAAAA
  
```

Restriction Sites: RsrII-NotI


[View online »](#)

ACCN:	NM_025942
Insert Size:	1191 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC011318</u> , <u>AAH11318</u>
RefSeq Size:	1692 bp
RefSeq ORF:	1191 bp
Locus ID:	67059
UniProt ID:	<u>Q9CZ30</u>
Cytogenetics:	2 C3
Gene Summary:	Hydrolyzes ATP, and can also hydrolyze GTP with lower efficiency. Has lower affinity for GTP. [UniProtKB/Swiss-Prot Function]