

Product datasheet for **MC203808**

Nat8f1 (NM_023160) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nat8f1 (NM_023160) Mouse Untagged Clone
Tag: Tag Free
Symbol: Nat8f1
Synonyms: 1110002I11Rik; A1266962; Cml1
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC010796
 CCCACGCGTCCGCCACGCGTCCGCAAGGAGAAGTAACTCTAATTTGTCGGGCTTCGGAGGTGGAAAAG
 CCCCCACTGGTCGGGCCCTAGAAGCTGAGGGTTCAAGGAAGGTGTGCAAGGCAGGTATAGCTGTCTCTCCT
 GGATGCCAAGATTTGAGACCCAGAAGTCTCCCATGGTTCCCTATCACATCCGACAGTACCAGGACAGCGA
 CCATAAAAAGAGTCGTGGATGTGTTACCAAGGGCATGGAGGAGTACATTCCTCTACCTTTCGGCACATG
 CTTATGCTGCCCGAACCTCCTGCTCTTACTTGGGGTGGCCCTTGCCCTGGTCTGGTGTCTGGCTCCT
 GGATCCTGGCTGTTATTTGCATCTTCTTCTGCTCCTACTTCTGCGGCTCCTTGCCAGACAGCCCTGGAA
 GGAATATGTGGCCAAATGTTTGCAGACAGACATGGTTGACATACCAAGTCTTACCTGAATGTACATGGC
 GCCTGCTTCTGGGTGGCTGAGTCTGGGGGCAGGTGGTGGGCATAGTGGCTGCTCAGCCAGTCAAGGATC
 CTCCTACTAGGGAGGAAGCAGCTGCAGCTCTTTCGCTGTCTGTCTCCTCACAGCATCGAGGACAGGGGAT
 AGCGAAAGCGCTGACCAGAACTGCCTCCAGTTTGCAAGGGACCAGAGTTACAGTGATGTTGTCCTTGAG
 ACCAGCGCTTGCAGCAAGGTGCTGTGACTCTCTACCTGGGCATGGGCTTCAAGAAGACAGCCAGTACT
 TCATGAGTATATTCTGGAGTTAGCAGGATTTGTACAATTCAATTAAGTACTCCTTCCCTTCTGCCTA
 GGAGGGGTGGCTGTGACCTTATGCTCCTGTGCAGCAAGCACACTTCTCTGCACTCTGCTACAGGAACCAG
 TGAACCTGTGATGTCAGTGTGATTAACAATAAAAGTTTGTGGTGCACACCAGAAAAAAAAAAAAAAAAAAAA
 AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: RsrII-NotI
ACCN: NM_023160
Insert Size: 669 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC010796</u> , <u>AAH10796</u>
RefSeq Size:	1015 bp
RefSeq ORF:	669 bp
Locus ID:	66116
UniProt ID:	<u>Q9JIZ0</u>
Cytogenetics:	6 C3
Gene Summary:	May play a role in regulation of gastrulation.[UniProtKB/Swiss-Prot Function]