

Product datasheet for MC203827

Taz (NM_181516) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Taz (NM_181516) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Taz
Synonyms:	5031411C02Rik; 9130012G04Rik; AW107266; AW552613; G4.5
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC015305

```

CGACCCACGCGTCCGGTCAAGACTTGGTTTCTAGTGCGCATTGCTTTGTACAGCAGTCCCCTCTCGTCGC
CTGAGCAGCGGCTCTTCTCTATCCAGAGAACAGTCTAGGGAGGTGCTGAGCTGTTCCACCAGTGACCT
TAGTCTAGTCCTCTGCCAGTCTGTTGGAGATCCTTCTCCTGAGTGGTAGAGGACATCGGCTGGCCTGGGG
CAAGAGGAACTCTGGGCTCCGGGGGGGGGGCCGGAGCCGGAGGAGATGCCCTCCATGTGAAGTGGCCAT
TCCCCGCGGTGCCCGGCTCACCTGGACTCTAGCCAGCAGCGTCGTCATGGGCCTAGTTGGCACCTACAG
CTGCTTCTGGACCAAGTACATGAACCACCTCACCGTGCACAACAAGGAAGTGTGTATGAGCTCATTGAG
AACCGAGGCCCTGCCACCCCTTATCACCGTCTCCAACCACAGTCTTGCATGGATGACCTCATCTCT
GGGGGATCCTAAACTCCGCCACATCTGGAACCTGAAGTTGATGCGTTGGACCCAGCAGCTGCAGACAT
CTGCTTACCAAGGAGCTGCACTCCATTTCTCAGCTTGGGCAAATGTGTGCCTGTGTGTCGAGGAGAT
GGTGTCTACCAGAAAGGATGGACTTCATTTTGGAGAAGCTTAACCATGGGGACTGGGTGCACATCTTCC
CAGAAAGGAAAGTAAACATGAGTTCTGAGTTCTGCGCTTCAAATGGGGAATTGGACGGCTGATTGCTGA
GTGTCATCTGAATCCCATCATCCTGCCACTGTGGCATGTTGGAATGAATGATGTCCTCCCTAACAGCCCA
CCCTACTTCCCCCGCTTTGGACAGAAAATCACCGTGTGATTGGGAAGCCCTTACGTACACTCCCTGTGC
TTGAGCGGCTCCGGGCAGAGAACAAGTCAGCTGTGGAGATGCGGAAAGCCCTCACAGACTTCATTCAAGA
GGAATTCAGCGGCTGAAGATGCAGGCAGAACAGCTCCACAATCATTTCCAGCCTGGAAGATAGGCTTTG
CCTGCCAGTCTGTGCCCTGCAGTTCTTAGCCAAAGGAACATTGGGTCTTCAGGATTTGGGAAGTGGT
GTACAGGAACCTGACCAAGCTTGGCTCTTTACTGCTTTTGGAGATTCTCCCTGTACAAAACCTGAGG
CTGGGAAATGGACTGATGCTTTTAAATTCAGATAGATTGGCTCCCAATTCCTGCTGGATACCCTTTTTAT
TCTTGAACATCTGACTCCTCTGGGTCTTGAGCTAAAAGAAGAGCCTTGGCTGTCCCTTATACTTGGTGGA
GAAGGAAGGGACTCTTAGGCAGGGGGCCCTTCTCACCTCTCAGGGAGGGAACAGAGCAGGTAGTCAA
GCAAAATAGTTGGGTAGGAGAACAGGTTTCTTGGTCATAGGTACAGAAAGAAGACTCAATTGGTCCCTTGG
ATGGGCTACTACTTGCCCACTGCTCACAAACCTCATCAGGGATAGCTGTGAGCCAGGTCATTCTACCTGA
GCTACTGTACAAGTTTGTCTTGGGTGGCCTGTAGACAAGAAGAGAAGGTATTGGCCAGGATTCAAGC
ACAGCTAAACTTAACTAGATGTTTCTTGTATATAGCAGAAGAAGGACCCCAAGACCTGACCATTAACCC
CCAGTGCAGGGCTACCTGCTTCTAATCTTTTTTGCCTTGTATCCTCTATATCCTGGGATCTGGGT
TTAATGTTTTTAGAGTAATAAAACATTTGTGTTGTGCCACCAAAAAAAAAAAAAAAAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_181516
Insert Size:	789 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC015305 , AAH15305
RefSeq Size:	1807 bp
RefSeq ORF:	789 bp
Locus ID:	66826
UniProt ID:	Q91WF0
Cytogenetics:	X 37.95 cM
Gene Summary:	<p>This gene encodes a mitochondrial phospholipid-lysophospholipid transacylase necessary for normal composition and content of cardiolipin. In humans, mutations of this gene result in Barth syndrome, most often characterized by cardioskeletal myopathy, neutropenia and abnormal mitochondria. This gene is distinct from the gene encoding transcriptional coactivator with PDZ binding motif. Both genes share the gene symbol Taz. Multiple transcript variants encoding different isoforms have been described. [provided by RefSeq, Mar 2010]</p> <p>Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 3' coding region, compared to variant 1. This results in a shorter protein (isoform 2) compared to isoform 1. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>