

Product datasheet for **MC203850**

Ginm1 (NM_145418) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ginm1 (NM_145418) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ginm1
Synonyms:	AA589616
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC013529 CCCACGCGTCCGCTCCC GCCCGCCGGCCGGTGTGTCCAAGATGGAGGGCCGCTGTCCGGCCCGCTG TCCGTCGGTTGTTGCTGTTTCATAGCGCTGCCAGCCGCGGGATGGCTGACCACCAACGCGCCAGGCCG CGTCCACAGCCCCGAGAAATGGCATCCAAATCAATGTAACCTACCCTGAGCAAGAGTGGGAAGAGTCTGA AGAACAGTTGTTCTTAACATAACCTATGAGCGTGGACAGGTATATGTAATGACTTACCTGTAAATAGT GGTGTAACCCGAATAAGCTGTCAGACTTTGATAGTGAAGAGTGAAAACTGGAAAAATTGGAGGAGAAAC ACTATTTTGAATTGTCAGTGTGAGGATCTTAGTTCTCGAGAGGCCTGTGACGTACAGTGCCAGCTCCCA GCTGATTGTCATCCAAGGAGAGTTGTGGAGATTGACGGGAGACAAGCTCAACAAAAGAATGTTACTGAA ATTGACATTTTAGTTAAGAACCAGAGAGTACTCAGATATTCAAGCTATTTCTTCTTTGGAAGAAAGCA TGCTTTATTCTATTTCTCAAGACAGTGACATTTTATTCACGCTTCTGACTTTTCCAAAAAGGAAGTGT TAGCTCACTGCAGACCACGACTACCTCATGGGGAATGTGGAAACCACTGTGGATGGAAATGCTTTA CCTGGCAAGTTGCCTGAGACTCCTTTAGAGCAGAGCCTCCATCTTCATACAAGTAATGTGCCAGTGGA TGGAGAAGCTGAGAAAAGCGCTCTGTAGTTCTGGAGCAGTGTGGTCCCTGTGCTTTCATGTTTCTGGA CGTCATGGTGGTTGGAGTTCTGGGAGCAGCAGGAGTATAGCAGTGTAAAGCTGCTTCCCAGTGTGT GAGAACAAGGAATTCTCAGGTGGATAAAATGAATGGGATCTCGGTACCTATCATCTTATACCCAGACG GCTCTGAGAAAACAGCACAAAACCTTACAGATAAAACAGATATTTAAAGCACCATCTCCAATCATGGACT CTATTGCCCCAGACTTGCCATTTGAATACCATTTTATTTAAATCATTCAACTTAAACAAAAGAAATTGCC AAATCCCAAGACCTTTACAGTTGGCCCTTTAGTGTGAGGCTCCGTGCACAGGCTGGCTGTTTCAAGA TGAATCACTAAGGAAAGTTTAAATGAAAAGAAGATGACTGAGAAATAACCAGACTGAAAGGGAAGCCAA ACCTCCAGGTGCTTTGCTGCAGAACTGTATCTATTTGATTACTATTGTAAGATTTTATAAGAAGCAC TTTATGAAATTCTAATTTAAAAGCTGGAGAAGCTAACATTTATTTTTTTTTTTACTCAGTACAATGTGCTT TTATATCACATACATGACACAGTATAAAAATCACTGAAGTGTATCATGATTTGTATCTGTACCAAAACCC CTTACCTTGAATGCACCAGGATTTATAATAGTGACATGGTCTCACTGTAAAAATAACTTTTTTAAAAAAA AAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_145418



[View online »](#)

Insert Size:	984 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC013529</u> , <u>AAH13529</u>
RefSeq Size:	1548 bp
RefSeq ORF:	984 bp
Locus ID:	215751
UniProt ID:	<u>Q91WR6</u>
Cytogenetics:	10 A1