

Product datasheet for **MC203858**

Sectm1b (NM_026907) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sectm1b (NM_026907) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Sectm1b
Synonyms:	1810003C24Rik; K12; Sectm1
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC010805 GCAGGTAAGAAGCTGGAGCACACATGCAGTGATTTGTCCAGAGTAAAGGTGAGGCTGAGGAGTGCCAGG CACCTGGACATCAGGAGCAGAGTGGATCTGTGTCTCCAGGTGCTCCACCAAGCTCCAGGAGGAGC AGGTGCAGCAGCTGGGTTCTTCGGTGTGAAGTGAATGGGAAGAAGAGAACAGAGCATTCTTCTGTGCAT CTGATCATGGTGTCCAACAACAACAACAACAACAACAACAACAACAAGACACTGACCAGAACAATTGT CTTCTGAGGCCATCCTAGGAATGAGCCATGCTGGCCTACTCTGTAAACATCCTCTGGCCTGTTTCCAGA ATGCTCTGGGCCCTCCTTCTACTGGCGCCTCCCTGAATGCCATAACGATGTCTGGGACGAACCTTGTT GCACTGAGCATGAAGTATCTGTAACAGAGGCAGCCGTGTGGTATGGCCTGTAATATCTCCAACAATCT CAGAGACGTCACCATTGAGTTGGTTACCAGTGAAGAGACTAGCATCATCTTCAATCATACGCTCCAGGA AACTACTTAAGGATTCATGGCAGCTTCATATCAAGGAGTCCAGGCCAGCTGGTGATCAGATGCTC AGGACAAACACTCAGGGAACACTCCTGGAAGCTGCATGGATTCCAGGCAGAGTTCAAAAACCTTCAACCT GACTGTTAATGCCGACAGACAGAGAAGACAGAGGACTTGCCAGTCACTAAGGTCCTGATAAGCCCCCA ACTGCAGTGCCGACAGAAGTTATCATCATCATCGCCATTGCTACCACCATCATCATCACAGGAATCGGTG TATTTCGTTTGGTACAAGCAATCCCTGTAGCTCCACAGATACAAATGTCAGTACCTTGCTGATTGATG GTCTCCTGGTATCCCCTACCTGACACTGCCCCCTAAAGTCTCCCAATCCTGCCTCAAAAGATACTCAA GGGACCTTAGTAGCTCCAACATCAGGCTACAGAGAACCTCACCTGGCTATGGCTGACAAATCAAGCCCT GATTGACAAGATCACAGGAACCTCTGTCTCCCAAGCCCTCATTTGACCTTGGGACAGGTAGACGTCAG AATGGGGTCTGTACTCCTCAGGTTCTCCCAAGACTCCAGGACAATTTGCACTGAGGACGGCACTGTGGT TTTTAACTGAGGATACCCTGGAAGTTGTCTGAACCCCTCCTGTACATCAAAGTATTTTCAGAGTATCCCAA AAAAAAAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_026907
Insert Size:	639 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC010805</u> , <u>AAH10805</u>
RefSeq Size:	1273 bp
RefSeq ORF:	639 bp
Locus ID:	58210
Cytogenetics:	11 E2