

Product datasheet for MC203920

Cpsf4l (NM_026682) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cpsf4l (NM_026682) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cpsf4l
Synonyms:	0610010C04Rik; 1500000C01Rik; D11Ertd636e
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC016190

```

CTCTTTCTGCTTTTCCGGCGGGAGACGGAAGATGGAGGAGTTCATTGCCGGGCTGCAAGGGTTCACGTTC
GCCTTTGAGCTGGATGTGGAATCGCAGAAGGGCACTGGGCTCCTGCCTTTCCAGGGCATGGACAAGTCCA
GTTTCAGCTGTGTGCAACTTCTTCGCTAAAGGGCTCTGTGTGAAAGGGATGCTATGCCCACTCCGGCACGA
GCAGGGTGAGAAGCTGGTGGTGTGCAAACTGGCTTCGAGGCCTCTGCCGGAAGAGTGATTGCTGCGAC
TTCCTGCACCAAGTATGACGTGAGCAAGATGCCCGTGTGCTATTTCCACTCCAAGTTCGGGAAGTGCAGCA
ACAAGGAATGCCTGTTCTGCACTTGAAGCCAGTCTCAAGCTCCAGGACTGCCCGTGGTATAACCAAGG
GTTCTGCAAGGAAGTAGGTCCCCTGTGTAATACCGCCATGTTTCATCAAGTACTGTGTCCCACTACTTC
ACCGGCTTCTGCCCCGAGGGACCTCAGTGCCAATTTGGGCACCCTAAGATGAGTCCTCCGTTCCATCCCA
GTAATGTGAAGCTGCAGCCTGTCAACCAGCCATGGGATTCAACAATTCTACTGGCAACTCTGGTGTC
CCCTTTCCAGGCAAAGCCCATGGTCCACGGTCAGAGGAGGTGGAGTCTGCCTCAGGCCTGCAGCTCAAGG
GCACACCTGCACCTTAGGCCACCATGCCAGCAAAGTAAGCTGTACCCGTTGGAGCCATCAAAACCTAGT
CTGAAGCGAGGGAGTTGGACGTTCTCCTAGGTAAGGGCTGGCTCCATGGCCTCCTCAGGGTATTCTGTAC
CCAGTGGTCCAGTGAGCTCCATCACGGTAATGAGATGCAGGATGCAGAAAGAGCCTCCTCTACCTGGCTT
ATTTCCATTTATCCTGCACATGCCCGAGATGCCAATCATACAGCAGGATACGAGGTGAGTGTGGA
AGGGAGGAGTCATTTCATGAGCCAGGAGAACAAAGTGCACCTTGAAGCTGTGAAGTCTACAGAGTCAGTCAG
TCAACATCCCATGCTTGGCTCCCACGCTCCCGCAGACTGCTATGGGTCAACGGCCTCAAGTGCTCCCAAA
AACGGGCAATCCGATCCAACGATTGCTTTTACTTTGACAGCTGTTGCCTGTACAGCTTATTATACTCG
AGAATATACTGTGTGTTCTGTGTAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
  
```

Restriction Sites:	RsrII-NotI
ACCN:	NM_026682
Insert Size:	531 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC016190 , AAH16190
RefSeq Size:	1240 bp
RefSeq ORF:	531 bp
Locus ID:	52670
Cytogenetics:	11 79.17 cM