

Product datasheet for MC203932

C030006K1IRik (NM_176828) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	C030006K1IRik
Synonyms:	11100331I3Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC019462

```

GAAGGATGTGGCCGCCTAGTGGGGCAGTGCCTAACTTGGCTCTAGTCTTGGCGAGGTCGCAGAGAGCTCG
GACCTGCAGTGGGGTCGAGCGTGTCTCCTACACACAGGGTCAAAGCCCGAGCCAGAACCCGCGAGTAT
TTCTACTATGTGGATCACCAAGGCCAGCTTTTCTAGATGACTCCAAATGAAGAACTTTATCACCTGCT
TCAAAGACCTGCAGTTCTGTGCTCTTCTCCCGTCTGAGACCAACCACAGCGGGCGTACGAGGC
TTCCTTCCCTTCTCTCGCTTTCGGCAGAGAGCGTAACTTCTGCGCTGTGAAGACCGGCCAGTGGT
TTCACGCACCTCCTGGCTTCCGACAGCGAGTCACCTCGTCTCTCCTACTGTGGCGGTGGCGAGGCGCTGG
CGATACCTTCGAGCCGGCTCGCTTGTACCCCTGGCGGCCAATGGGCGCTGTACACCCGGCCCCGGA
GGCGCAGGCGGTGTGGGTCTAGTGGCTCGGCCCTGGCTTCGAGCTCAGCGCTTGCTTCGAGTATGGG
CCTAGCTCCCTACGGTCCCTCTCACGTACACTGGCAAGTCCGCGCATCGCCCTAACTATGGATCTGG
CCCCGTGCTGCCGCTGCCCGCCGCCCTGAGACTCCAGCTGGGGCTGTGTGGGAGGGTTAGAGCAAA
CTCTGGGTTTTGTAGTACTCGCAGGTTACCCCTCTGCTCCAACACTCTAGATTAACTTAGATGTCCGCTG
TCACCTGCACGCTGGTACTGAGAGCGCTGCGCGCAGAGGTGAGCCAGGGTGGTTGGCGCTGTCCAAG
GTTCTGGTGTGCCCTGCTATCCTCTGTCCCTACAGAGGTCTCTCTCGCCACCTTAACCCCTTACACC
GGCTTTGCGTCCACGGATTGGCTGCGCTGGTCAGTCAGGCGAGGCGGAGGCGCAAGGCGCGCTGG
GCCGCTGCGCGTGCAGCCACCCCAACCCCCCTCCGCTTCTCAGTCTTGACGCGCAAGCACCGC
GGGGAGTGCCTCGCCACTGGAGCAGGTGAGGGGTGGGGGAGATCCTGGCCCAAGGTACCCGGGCTTG
GGGCGGGATCGGACCCCTCTAGGTTCTGACTGCCTACCCGAGAAAAGGTATAGGATCTGGCAGTGGAA
GGGATGGAAACACACCCTCAATGAGCAGGCTCGGACCCACTGGAACGGGAGGGGGCGTTAGGGGAGG
TAGTGTGCTTAAAGCGCCGATCGCAGGGGTGCCAGAGACCTCCGCTGGGCGGACAGAGGACACAG
GACTCTCCCAACCCTTCTCCGCTAGATAGGCATCCAAGAGACCGACCCCACTCTGCTCAGGATT
AGAGGTTCTGCTGATAAAAGGATCTGCTGGGCTCATCTAATAAGGGGCCCTTATCTGACATGG
CTGTCTGCTGCTGTCACATCTGGGCTCTGCATAAGCAGCAGCTCTGATCCCTAATGCACAGATCAAG
GCTTTGCTACTGAAGGCTGAACAGCTAGCTGTGCAAGGCTTGAGCCAGGACCAACCAAGTCATTCTGT
TCCAGGAAATAAACCCAGAACCTAGCAAGCGTCCCACTTCTTGGGCTCCACATGTCTGAGATCTA
GAACTGCCTTTCCTGCGCCAAAGCTATTCATTAATACTTACTGCTGTGCTGCCTCAAAAAAAAAA AAAAA

```



Restriction Sites:	RsrII-NotI
ACCN:	NM_176828
Insert Size:	657 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC019462</u> , <u>AAH19462</u>
RefSeq Size:	1755 bp
RefSeq ORF:	657 bp
Locus ID:	223665
UniProt ID:	<u>Q8VE95</u>
Cytogenetics:	15 D3