

Product datasheet for MC203976

Nipsnap3b (NM_025623) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nipsnap3b (NM_025623) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Nipsnap3b
Synonyms:	2700063N13Rik; Nipsnap3a
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC019505

```

CGTATCCGAGCGCCCGGCGGCGGTCTTCTGCGCGCTCAGGACCGCACCATGCTCGCGCTCCGAAGCGGC
CTAAGGACAGCGCTGGCCCCGCGGTTCTGACGCCTCAGGTATGTTACCTTTTGCTACAGGCCCCAGGC
AAAGCAATGGCACATTCTATGAATTCGACCTATTTTCTCAAGCCCTCAAAGACAAATGAGTTCCTGGA
AAATTTTAAAGAACAGTGTTCACCTTCGGACAGCTCACTCTGAAATGATTGGCTATTGGACTGTAGATTT
GGAGGCAGAACGAACAGAGTTTTCCATATTTGGAAGTATGATAATTTTGCTCATCGAACTGCAGTCCGCA
AAGCCTTGGCCAAGGATAAGGAATGGCAGGAACGATTCTCATTCCAAATTTGGCTTTCATTGATAAACA
AGAGGTTGAGATTACCTACTTGGTACCATGGTGCAAAATAGGAACGCCTCCCAAGGAAGGAGTCTATGAA
CTGGCTACCTTTCAGATGAAGCCTGGCGGCCAGCTCTGTGGGTAATGCATTTAAAGGGCAGTAAATG
CCCATGTGGAGCTTGGCTACTCTACACTAGTTGGTGTTCACACCGAATATGGAGCCCTCAACAGAGT
TCATGTTCTGTGGTGGAACGAGAGTGCAGACAGCCGAGCAGCTGGGAGACACTGGTCTCACGAGGACCCT
AGAGTGGTGGCTGCTGTTCTGGGAGAGTGTCACTTACCTAGAGTCCCAGCAGAATACGTTCTCTGATTCCAA
CATCATTTTACCACCTGAAGTAGTTTCTACCCAGCCCAGAGCATCCCATCACAGAGGTGAGCTGTATCC
ACAGCGGGCACTTCTGTCTCCAAGATGGAGTCGCCCAGTTTATTGTCTGAAGGACGGGCTGAGCTATTA
CACTGTGACCTCTGCATCTTGAAGGCTACCTGAGTCATTTACTACGAGTAGTTCAGGAAATGTTTCTGTT
CATTTTCATAACTACTGGATATCAGAAAGCATTGTTACTCCGTGCCCTGGTACTTCATTGTTTACAATA
TCTTGCCCTTCTGAGATTTTTATGACTAAGTCTTGATTTCAACATGATTTTGAATTTTTTCCCTCAAAG
TATCTGTAATTAATAAATAACAATGTCTGCCTTTTTATCATGTTTCAATTTTCTGAAATTTAAACTAT
TTTTTACCAGTAATCCTTTCTTTGTGATAATAAGTGGGATTTATAGAAATATTAATCTAGTATTTGTCT
TTATACTTTAGAAAAGTATTAATTATGGTTTCACATTCATTATTAAGTCAAAATTGGAAAATAATATATC
CTTGAGACTTGGTGATAATACGGTCTTAGTGTGTTTGTGTTTAAATTGCTGTACTACTTTTAACTTTAC
CTTGAATTACTCATAAATATTATGAGACACTTATTTATTTCAATGAGAAAATGGAATCTCTGACCTAAT
AAAAGCCTTCACATTCATTCTGAAAAAAAAAAAAA
  
```

Restriction Sites:	RsrII-NotI
ACCN:	NM_025623


[View online »](#)

Insert Size:	744 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC019505 , AAH19505
RefSeq Size:	1507 bp
RefSeq ORF:	744 bp
Locus ID:	66536
UniProt ID:	Q9CQE1
Cytogenetics:	4 B2