

Product datasheet for **MC203986**

Slc25a35 (NM_028048) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc25a35 (NM_028048) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Slc25a35
Synonyms:	1810012H11Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC019996

```
GCTCTAGCTCCCGCAGCTACTTTGCCACGCTTATTCTTAGCCGAATCCTAACTTCTTTCACAACCCCTCC
TCTCGCTGTTGTGCGAGCCACCTTCGGCTTTTCCCTCTGATTCTTTTCTGTTGTTGCTTTTCTTTTCCCT
GCTTCTCCAGCTACCATTTTCGTAACCTCCGTTTTCTAAGATACGCTCTGGTACCCCATCCTAGGCATTC
CCACCTGCCAGCTTCTAACTTCTGCTCCTCCCTGATCTCAGTTCCCTTGTCAACAGGGCCTTGTTGCCAT
ACCCTACTACACTCCTCTTTCTTTCTTTTCTGCTACCTTCTGATCCCCAACTCCCTCTAACCTCACTTGA
CTTCTTGACATTTGACCCAGACAATTTAGCTCTCTTAAATTTCTACCTGTTGCAGCCTGGGCTTCAAAGA
GCCAGGTGTTGTTTCATATCCATCCCTGCTTTCTGAAGCTGTTTCTACCTTTACCATGGACTTCTTAATGA
GTGGTGTAGCAGCCTGTGGAGCCTGTGTGTTACCAATCCCTTGGAGGTTGTGAAGACCAGAATGCAATT
GCAGGGAGAACTGCAGGGCCAGGCACCTACCAGCGCCACTACCGAAACGTCTTCCACGCCTTCTTCACC
ATTGGCAAGGTGGATGGCTAGCTGCCCTGCAGAAAGGGCCTGGGCCCTGCGCTCTTGTACCACTTCTTGA
TGAATGGCATCCGACTGGGCACCTACGGGCTGGCAGAGTCTCGGGGTTACCTGCACACTAATGAAGGCAC
CCACAGTCTGTCCGCAGCGCCGAGCTGGGGCCCTTGCTGGGGTTCATGGGAGCATATCTGGGAAGCCCA
ATCTACATGGTGAAGACACACCTACAGGCACAGGCAGCCTCTGAAATTGCTGTAGGGCACCAGTATAAGC
ATCAGGGCATGTTTCAGGCACCTAACCGAGATTGGCCAGAAACATGGTCTAGTGGGGTTGTGGCGTGGGGC
CGTGGGCGGCCTGCCAGAGTTGTATCGGCTCCTCTACCCAGCTGTGCACCTTCTCATCTATCAAGGAC
CTCTTGAGCCAGTGGGAGATCTTCCCTCCCAGAGCTGGAAGGTGGCTCTGGCCGCTGCCATGGTGAGTG
GTGTCGCCATAGTCGTGGCAATGACACCTTTGATGTAGCCAGCACAAGACTTTATAACCAACCCACTGA
TACTCGAGGCAAGGGCCTCATGTACCGGGGATCTTGGACGCTCTGCTGCAGACAGCGCGGACAGAGGGC
TTTTTTGGCATGTACAAGGGTATAGGTGCCTCCTATTTTCGCCTCGGCCCCACACTATTCTCTCCCTCT
TCTTCTGGGATCAGCTGCGCTCCTTTTACAACACATATGCTAAATAACAATCACTCTCTTGTAACTCCAG
ATCAACACTCTTTGGGCACCTCTGCCTCTACTTCTATCACAGTACTGCTGACAGATTGACCTCTCAAT
GCCCCAAGGGGGACAGCCACGCCCACTCCTCTGTTCTCCTAGGGTTGAGTCACTGAGTGGGTCTTGCTT
TAAATTTTCCACAAATCCCCGATATATTTTATACCATCTCAGGGTTGAATCAGTCACTCCAAAGTAT
ATTTCTACTGCCATTCTTGGGTTCCCCCAACTCCCTGCTCTATACACAGCTTTAACTCCCCCCCCCTTC
CAAGACCAAAGGGAAGTGGGAGACCCAGGTGTCCTCTTAGATTTGAAGGCACAGAAATTATGTTCTTCAT
AAAGCAAGTTTATTTCTGCAAAAAAAAAAAAAAAAAA
```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_028048
Insert Size:	903 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>BC019996</u> , <u>AAH19996</u>
RefSeq Size:	1786 bp
RefSeq ORF:	903 bp
Locus ID:	71998
UniProt ID:	<u>Q5SWT3</u>
Cytogenetics:	11 B3