

## Product datasheet for **MC204009**

### Dnase111 (NM\_027109) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Dnase111 (NM\_027109) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Dnase111  
**Synonyms:** 2310005K03Rik; AI324885; Dnase11l; Dnase11; Dnase11; Dnase11; Dnase11; G4.8  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >BC023246

```
CGATAGAGCATTTCCTAGTATTTCATGTGGTACCAGAACAATAAATTTCTGCCATGAGGCATTTTTCTA
AACTGCCACTCCAGTGGGACCTGAACAGGTGGTCAAGCTCTGTGTACAGAAGTTGAAATACGTTCCAGG
AAGAGATTCAAGAAATGTTCCAGGGCCAACAGCAGGGAGAATACTAATTAGCTGAGATCAGCTGAGTGCC
AATGCCCTTTGGACAACCTGGATTCTTTGGAGGGTTCCTGATGCACACATAGCAATGCGTGGCCTTGTT
ATGGCTCCCTGCTCATTCTCCTGGTGGTGGGACTGAAGCCTTCGTATCTGTGCCTTCAATGCCACA
GACTGACACTGGCCAAGCTAACCAAGGAGTCAGTGATGGATACCTTAGTTCAGATCCTAGCCGATGTGA
TATCATGGTGCTTCAGGAGGTAGTAGATTCTTCCAGAATACTGTCCCCTTTCTGCTTAAAAAATAAA
AGTTCCAGGAGTTACAGTTCCTAACAGCTCATTGCTGGGGCGCAGCACATAACAAGAAAAGTATGTGT
ATATCTACAGGTCTGACAAGACACAGGTCTAAATTTCTACCAATACAACGACACAGATGATATTTTTGC
CCGAGAACCATTTGTGGCCATTTCACTCTCCCTAGCAAACTCTCCAAGTGTGGTGTCTGGTTCCTCTC
CACACCCTCCCAAGGATGTTGAGAAGGAGCTGAATGCCCTCTATGATGTGTTTCTGGATGTCTACCAGC
GCTGGCAGAATGAGAATGTGATTCTGCTTGGAGACTTCAATGCAGACTGTGCATCGCTGACAAAAAACG
TCTGAAAAGCCTGCTACTCCGGACTAAGGCAGGCTTCCACTGGGTGATTCTGATGGGGAGGATACTACA
GTGCGGGCCAGTACTAAGTGTACCTATGACCGAATTGTGGTGCATGGACAGGGCTGCCAAATGCTGCTGA
AGGCTGCAGCTACCTTTGACTTTCCCAAGAGTTTCAGCTGACTGAGGAAGAGGCTCTCCGCATCAGTGA
CCATTATCCTGTAGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGAAGTGA
CAGATGTGAAAAGTCAATGACAGCTTTTTTAAATAGCTGCTTGTCTATTCTGAGAAAAGTGAAGTGAAGT
ATAACTCAGCAACAAGTGCCTTGGGCTTGAATCCTGTTGGGTAGACAGGACTTCACTGAGGTAGTTAGG
AACAAAAGTCCAGTAGAAATGGGCTCAGATCATACCTAAGATTTCTGCAACTGATACAAATGAGCTCATA
CTGAAAAACCTTAGGTCTAAGTAATCTTGCTTCTCCATGCACTATACTACATGGCTCTGCCAGTGTGTG
CAGGTTTTAAATATTCAATAAAATCATTTTAACCTTTAAAAAATAAATAAATAAATAAATAAATAAATAA
```

**Restriction Sites:** RsrII-NotI  
**ACCN:** NM\_027109



[View online »](#)

<b>Insert Size:</b>	945 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC023246</a></u> , <u><a href="#">AAH23246</a></u>
<b>RefSeq Size:</b>	1535 bp
<b>RefSeq ORF:</b>	945 bp
<b>Locus ID:</b>	69537
<b>UniProt ID:</b>	<u><a href="#">Q9D7J6</a></u>
<b>Cytogenetics:</b>	X 37.94 cM