

## Product datasheet for **MC204126**

### **MIIt11 (NM\_019914) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** MIIt11 (NM\_019914) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** MIIt11  
**Synonyms:** Af1q; AI839562; Zfp692  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >BC022963

```

GTCTTCGGAACACCCCTAGTCATTTCATTGCTAACAGGGGATCATGAGGGACCCTGTAAGTAGCCAGTACAG
CTCCTTTCTTTTCTGGAGGATGCCCATCCCAGAAGTGGATCTGTGCGAGCTGGAAGGCCCTGGCCCTGTCA
GATACACCTACCTACGAGAGCAAGGATAGCAGCAGCGTTGGCAAAATGAACGGGCAGGCATCTGGGACGG
AGCAGAAGAACCCTGAAGGTGACCCCTCCTTGAGTACAGCACCTTCAACTTCTGGAGAGCGCCTATTGC
CAGCATCCACTCTGTCGACCTGGACTTGCTTTAAGGCCAAGATCTCTCTCCACCCCGTACCTTCATTAT
CTTCCCTCTAAGTCCCTCCTCTCTAAGCCTTCCACTTAACTTGTCTCTTCCCTCTATTGTGCTGG
TGGTGCTGATGAATCTGCCAATTGAATTGTAGGTTTGTATGTATGTACGTACTTATTTATCTATCTATAT
TTATTTATTTATTTATTTATCTGCCTACTCCGTTTCTCTCAAAGCCCTTAGGCACAAAGTCAAGGGTT
TAATCAATGGAGTACTAGATCATAGGATTTCTCCCGTCAGCCCGGAGAATATTAATAATCCCAGAAAA
GAATCACCTTATTTGAGGCACCACAGAGAGGTACAGTTGTGCAAAATGGGAGACCAATTATTTTGTGTTG
TTTGTTACATTTAATCAACTTCAAAGATTACTTGAGCCTCCCTAATGTCTCTACCAGGCCAGTAAAGCA
AACTACTGCTTCAAGAAATGCTGCTGCTGCTTACAATGCATTCTCTTCCCAGAAAACACGCAGGCAGTTAAAG
GTGGGAAACAAAGGAGTGGGATGAGGAAAAACATCGCACCGTTTCAGAGGTAGCTCCTTTAATATCCAGC
TCTGCTCCAGCTGTGGGTAGAAATGCTCAGAAAATCAGATTTGAGTCTTACAGGGAACTTCACTATG
AGGATAATTCTAAAAGACAAAAGATGTATCTGATCTAAAGAAACACCAGAGTATTATTAGGGGAATTTA
TATAAAATTGTGATTTTGTGATGGAATCAGTGATTTCTGACACAGTTCTTTGGCCGCCATAAGACCGTGAAGTA
AGCAGTTATGTTACAGACCCTACCCGTCCCATTCCTTCCAAAGCTATTGCAGATGAATTGCGGAGCTGTT
TATTCATTTATGGAATAAAGACGTTAAATCCTAAAAA

```

**Restriction Sites:** RsrII-NotI  
**ACCN:** NM\_019914  
**Insert Size:** 273 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC022963</a> , <a href="#">AAH22963</a>
<b>RefSeq Size:</b>	1392 bp
<b>RefSeq ORF:</b>	273 bp
<b>Locus ID:</b>	56772
<b>UniProt ID:</b>	<a href="#">P97783</a>
<b>Cytogenetics:</b>	3 F2.1
<b>Gene Summary:</b>	Cofactor for the transcription factor TCF7. Involved in regulation of lymphoid development by driving multipotent hematopoietic progenitor cells towards a T-cell fate.[UniProtKB/Swiss-Prot Function]