

Product datasheet for MC204171

1700011H14Rik (NM_025956) Mouse Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | 1700011H14Rik (NM_025956) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | 1700011H14Rik |
| Synonyms: | Ccdc198; RP24-142A23.1 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | PCMV6-Kan/Neo (PCMV6KN) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >BC026534 GGGTACAGAACAGGGAAGATGTCATGAAGGAAACACCTGCTTTGGAACATGGGCCTGGGCCACTCCAAGG CTCACCCCTCGGGTGATCAAAGTCACACCTTTACAAAGCCAAGAAGCAGAGACACCCTCAACTGGCCCTGT GTTCTTTGCATTGAATCGGAACCTGGAAGAAGAAAGCTCATTACAAGACTGCAGACCAGAGTCGAACT CGAGAAGGACAGCTACCACCGCTTCGAGAACTTGGTATGGAAGACTCCCGGCAGTATCCAGGGCCATGT ATCTTGACATCCCACTGAAACATGAAGAAACAAGTATTATCAAAAGGCATCCACCCGAAGAATTCAAAA GCTTGAACCCATTGACCTGCCACAAGCAATTACCTCTGAAAGGCTTCTGTGCCATCAAGAAGGCAGAACA AAGTCCAACACAAAGCAGGAAAATGAGAAGAAAATCAACTTCTATGTACACTTCTGGGAAAAGACAAT ATTTGCATAAGATGAAAATGCTAGAGATGAACCATAAAAGACAAGAGGCCCAAATGGAGTTGAGGAAGAG TCTTCTGAGTAAGGCGATGCTTGATATGCAAAAAGCTGAAGGACCATAATGGCAATAAAATCGCCCAGAGC AAGCCAAGGAGCAACAGCTATGACATTCTACCATACTGCCTAATGAAAACATAAACAGAGATCCAGGAA ATCCCCAAGATGAAGAATTCTTGGACTGTCACACAGAGAATGACTATTATGTCCGGAAAATTTGGCAAAAT GGAAACATGGCTCCGTGAACAAGAAGCCCGGGGACAGCTTTTCTGGGACAGTTCCAGCTCCGACTCTGAT GAGCTTGAGAAAGATGAGAGGAGACCACAAGCACTGGTGAGGACCAAGACAGAGAAGATCCCTCTTTATG ATGACTTTTATGGCAGTGCTTAAGAAAGCTAGACGCTACCCAGGTAGTAACGGCACAGAATGCCTTGGG CTTGTTTTCTCCATAAAATCATTATTTGCATATGGACTCCTCTGAGGGATAAAGCTATTTTGGTA TGCAGAAAGGATGTCATCTGGAATCTTGTGGTGAATAGCCTTGTATGTATGCCTTAAGTCTATGCCTTA AAACTAAGAGGCTGAAACAGAA A |
| Restriction Sites: | RsrII-NotI |
| ACCN: | NM_025956 |
| Insert Size: | 885 bp |



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>BC026534</u> , <u>AAH26534</u> |
| RefSeq Size: | 1191 bp |
| RefSeq ORF: | 885 bp |
| Locus ID: | 67082 |
| UniProt ID: | <u>Q9CPZ1</u> |
| Cytogenetics: | 14 C1 |