

Product datasheet for **MC204200**

Ces2a (NM_133960) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ces2a (NM_133960) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ces2a
Synonyms:	9130231C15Rik; AI266984; Ces6
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC025537

```

CCACGCGTCCGCACTCTCATGTGCAGCAGCACAGGTTTAGGCCTGAGATTCAACCATGCCATTGGCTAGA
CTTCCTGGTTGGCTGTGTGTTGTTGCCTGTGGGCTCCTGCTTCTCCTCCAGCATGTGCATGGCCAGGACT
CAGCCAGCCCCATCAGGAACACCCACAGGGTTCAGGTGAGAGGCAGCTTTGTCCAGTGAAGGACACTAA
ATCTGGTGTCCATGCCTTCTGGGAATCCCTTTGCCAAGCCTCCTGTAGGACTACTGCGCTTTGCACCT
CCTGAGGACCCTGAGCCATGGAGTGGTGTGAGAGATGGGACCTCACAGCCGGCCATGTGTCTGCAGCCTG
ATATAATGAATTTAGAGGATGCGAAGGAAATGAATCTGATCCTGCCTCCCATCTCTATGTCCGAGGACTG
CCTGTATCTCAACATCTACACACCAACCCACGCCAGGAGGGCTTAACCTGCCTGTGATGGTATGGATC
CACGGTGGTGGACTGGTGTAGGATCAGCTTCCATGAATGATGTATCTAACTGGCAGCCACTGAGGAAA
TAGTAATTGTGCGCATCCAGTATCGCTTGGGGTCCCTGGCTTCTCAGCACTGGAGACCAGCATGCCAG
AGGCAACTGGGGATACCTGGACCAAGTGGCTGCCTTACGCTGGGTGCAGAAGAATATTGCTTACTTTGGA
GGCAACCGTGACAGGGTCACTATATTTGGAGTCTCAGCAGGAGGTACAAGTGTGTCTTCTCATATTTCTGT
CCCCAATGTCAAAGGACTCTTTCATGGAGCCATTATGCAGAGTGGAGTGGCCCTGCTGCCTGACCTCAT
CTCTGACACTTCTGAGGTGGTCTACAAGACAGTAGCCAACCTATCTGGTTGTGAAGCCACAGACTCAGAA
GCCCTGATCCACTGCCTGAGAGCAAAGAGTAAACAAGAGATCCTGGCTATTAACCAGGTCTTCAAGATGA
TCCCTGCTGTGGTGGATGGGGAGTTCTACCAAGCATCCTCAGGAGCTGCTGACCTCTATGGATTTTCA
CCCTGTCCCAGCATCATTGGTGTCAACACAGATGAATGTGGCTGGGGAGTCCCCATGTTTCATGGGCCTT
GACCATATCATAAAGAACATAACCAGAGAGACCCTGCCAGTGTCTGAAGAACACAGCAGCACGGATGA
TGTTGCCTCCAGAATGTAGTACCTGCTAGTGGAGGAGTATATGGGGACACTGAGGACCCAGAGACCTT
GCAAGCACAGTTCAGAGAGATGCTGGGGACTTCATGTTTGTGATCCCTGCTCTCAAGTGGCACATTTTC
CAGCGTTCACAGCTCCTGTCTACTTCTATGAGTTCCAGCATCTGTCCAGCTTCAAGCATGTACAGGC
CATCACATGTGAAGGCAGACCATGGTGTGACGTTGCCTTTGTCTTTGGCTCCTATTTGTGGACATGAA
TTTGGACCTCACTGAGGAGGAGGAGCTGCTAAAGAGGATGATGATGAAGTACTGGGCAATTTTCGAAGA
AATGGGAACCCCAACAGTGAGGGTCTACCTCCTGGCCTGTGTTGGACATGATGAGCAGTACCTGCAGC
TGGACACCCAGCCTGCTGTGGTTCGAGCCCTGAAGGCCAGAAGGCTGCAGTTCTGGACCAAGACTCTGCC
CCAGAAGATCCAGGAGCTCAAGGGTCTCAGGACAAGCATGCAGAGCTGTAGCGTCTGTGCAGGGAAG
GTGTGTGAGATTAAGCAGATGGTGAATTTCTGAGGTTAAAAAAAAAAAAAAAAAAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_133960
Insert Size:	1677 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC025537</u> , <u>AAH25537</u>
RefSeq Size:	1805 bp
RefSeq ORF:	1677 bp
Locus ID:	102022
UniProt ID:	<u>Q8QZR3</u>
Cytogenetics:	8 D3
Gene Summary:	Carboxylesterases that catalyzes the hydrolysis of pyrethroids pesticides. Hydrolyzes permethrin faster than cypermethrin.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longer isoform (1).