

## Product datasheet for **MC204202**

### 1110059E24Rik (NM\_025423) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** 1110059E24Rik (NM\_025423) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** 1110059E24Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >BC023385

```
CTCAATCTGGGGCGGAACGTCTTTACTCTTTCTGTAACCTAAAGGCTGCACCGATTTACAGTGCTTGA  
GATTTTGGTGATGAGTTCCAGAAAGGCAACGTGACTCGGTCCCGCCTCAGAAGCACCAGAATACGTTCC  
ACCTTCAAAAATGACAAGTTTGATAAGAGTGTTACAGCAAGAAAATTAATGCAAACTTCATGATGGAG  
TATGTCAGCGCTGCAAAGAAGTCTTGAGTGGCGTGAAAGTACAGCAAATACAAGCCACTGTCAAAGCC  
TAAGAAGTGTGTTAAATGTTTACAAAAGACAGTAAAGGACTCTTACCACATAATGTGTAGACCATGTGCC  
TGTGAACCTTGAAGTTTGTGCCAAATGTGGAAAGAAAGAAGAGATCGTTATTCCGTTTAAATAAAGAACCAG  
ACGCTCCAGAAAACACTGAAAATGAAGTTCTGGCCATAGAAGGCGCTGCCGACGGAAGGAAGACAGTGA  
TGAGGACTTAGACGCTGAGTCTGACTCGGATGGCGAGGATGGAGACACTCAGGCCTGACAACACCAAGTC  
AATTGCAGAGGGAGTCAGAAGAGGCTCATGAGCAACATTGGTTATTTATAACATTTGAGTAACAGAATCT  
TAAAGGTTTTTACAGTAACACTTAAAGTGTATATACAGAATATGTGTGGGCACACCATGGCATGTTGTAT  
AAAACCTATAGAAATGTTTAGAAAAACATTATGTCAAATTATTAGTGTGAGAGTTGGTGACTGCATC  
AGCACTATCCATGAAAACTAAGATGACTAGCTCAGCGTACATGTCTGTAACTGACCTCGGGTATGTGA  
CTCGTTAGCAAGAGCCCCCTTACATCGCACAAAGCAAAACAGCTTATTTTATACACACTACACACAGT  
ACAATCATGTGTAATGATTTTCTTCAAAAATGTTTTATATGCTTCATTCCATTTTTTAACTTCTCC  
TAGTTTGGAAATATACACTTAGAAAAATGTGTAGACTGTATAGTTACTGTCTAGAAAGGAAAACCAACC  
TTCATACTGGAGTGGAGATACTTACAGTATCAGAAGCTCTGTGAAGTTCTCCCTCGAAAGGTTATCCTT  
CCCAGAATCCATCCTCCGGAATTATCACTTTATTGCCAGTGTTCACCTGTAACAGTGAAGTTAGT  
GTCACGTGATTATTGAAGTCTGTGCCTTGGGTTTAAAACCTAGAGCCACAGAATAGTGGCCAGGTGAT  
ATGATGGGCTTTTCCACACTCTAGCTAGAGTCATAATTTCAATGTTTACTGTGTTTTCTGTGTGCTTA  
ATCTCTACAAATGTATATATTACACACATGTGCATGTACATATGTATGTGTATATTTGCTGAGCCATTGA  
AATAAATGGGAAGCTGTGGTTTTGCAAAAAAAAAAAAAAAAAA
```

**Restriction Sites:** RsrII-NotI  
**ACCN:** NM\_025423  
**Insert Size:** 468 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC023385</a></u> , <u><a href="#">AAH23385</a></u>
<b>RefSeq Size:</b>	1440 bp
<b>RefSeq ORF:</b>	468 bp
<b>Locus ID:</b>	66206
<b>UniProt ID:</b>	<u><a href="#">Q9CQ90</a></u>
<b>Cytogenetics:</b>	19 B