

Product datasheet for MC204250

2610044O15Rik8 (NM_153780) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	2610044O15Rik8 (NM_153780) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	2610044O15Rik8
Synonyms:	2610044O15Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC023890

```
CCACGCGTCCGGCGGCTTGAATGGTTCTCTGTATCCTGCTTTGAGCGATGTTGTTGGTGCTTTGAGGG
GACCTCAAGTAAGCTCCGTTTTTGAACATGGATGCAGTGACTTATGATGATGTCATATGAACCTCACT
GAGGAAGAGTGGGATTTGCTGGATTCTTCCAGAAAAGACTCTATGAAGAGGTGATGTTGGAGACCTACC
AGAATCTCACTGATATAGGCTATAATTGGCAAGACCATCATATTGAAGAACATTGTCAAAGTTCTAGAAG
ACATGAAAGGCATGAAAGAAGTCATATTGGAGAGAAACCTATGAACGTAATCAATGTGGGGATTACAGT
TACATGATTGGATGACTCTCAGAAGAGACTTTGAACTTTGGAATTTAATATTTTGGAGCTGGTATAGA
TGATGGGGACTTTTAAAGTTGGACCTCATGATTTTTTCATTATGCTATGTTTAAAGTATGGTCCCATATA
TTCCTGTACTCCAACAAGCTTATGGGGCCAGGTAGTGAATGTGGGAGTTTAAATATGTTTGTACAGG
GAGTAGCTCTATTCAGAGGAGTGACCATGTTGGAGGAAGTGATCACTCTGGTGGTAGGCTTTGAAGCTC
TACCCAGTGCAGAACAGACCCTCTGGTAGGCCCTGAAAGACAGTCTCTCCTGGATGCCTTTGGATCA
AAATGCAGAACTCTCAACTCCTTGCTGCTGCTGATGTTGCCATGCTTCCTGCTATGATGAAAATGGACTA
AGCCTCTGAAACTGTAAGCCAGCCCAATTAAGGTTTGCTTTTATGAGAGTTGCCTTGGCCATGATGTC
TATTCCTAGCAATGGAATGCTAGGACAGAAACCTGTGAGTGTATTTAATATGGTGAAGCCTTTGCACA
ATTCTGTAGGCTTCATCATCATGAAAGAATATCTACTGGAGAGAAACCTTTTGAATGATTCACTATGGT
TAAGCTTCTTACTATGGTTAAGAAGGAATGACTTTGACAGATACAGTTATCATTCTTTCAATTTAGAAGG
TAAGCATCACACAGGAGTTTATGACAAGCTCATTCTTGGGCACTGTCTACATCTGTTAAATCAGAAGTTT
TGCTCTCAGCAATGTAATCTGTATCATGCAGGAAACAAGCCCAGGTCACCTTTCTGTATTAGTTCTCA
TGACTTTGGATTATGGAATTACAGAGTTTCTGCTCGTTAATTCAAAAGGAGTTGTCTAGCTTCTCTTT
AAGGAAAAGTGAAGTTTAAATAACTATATAAATAATAGAAATACGTATGTAATGTTATACAGCATTGCTA
CATTAATTTTATATGAAACATACACAACGTGATACATCAAAATGTTAAATGTATAAACCTGTTCTCCAAA
GCAAAAAAAAAAAAAA
```

Restriction Sites:	RsrII-NotI
ACCN:	NM_153780
Insert Size:	267 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC023890 , AAH23890
RefSeq Size:	1417 bp
RefSeq ORF:	267 bp
Locus ID:	72139
Cytogenetics:	8