

## Product datasheet for MC204479

### Taf7 (NM\_011901) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Taf7 (NM_011901) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Taf7
Synonyms:	55kDa; AI607308; BB007175; Taf2f; TAFII55
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC029673

```

GGAAAGTTGTCTCCGGTGTTCGGCTCTCTGAGAGGGAGCTTTCGCGCCGGGGTCTTAGCCGTCTCATT
GTTAGCGCATAGGTGCGCCTGCCGGCGGCCGAGAACGCCAAGAAACCGGAAGGGCGGACGCACCGCTCCC
TGCTGCTTTGCTCAAGCCGCGAAAAGCTAAGTTCACAGTCGCTGTTGAACTGAGCTCACAAAGGAATTA
GATGAGTAAGAACAAGACGATGCGCCTCATGAGCTAGAGAGCCAGTTTATCTTACGACTGCCCCAGAA
TATGCCGTACGGTGAGGAGGGCAGTGCAGTCTGGGCATGTCAACCTGAAAGACAAAAGTGCAGTTCAGT
TACACCCTGATGGGCGTCACGGCATTGTGAGAGTGGATCGAGTGCCTTGGCTGCAAAAATTGGTCGATTT
GCCTTGTGTTACGGAAAGCTTGAAGAACCATCGATAAGAAAACCTTTTACAAGACAGCTGATATCTCTCAG
ATGCTTGTAGCCACAGTAGATGGTGACCTATACCCTCCTGTGGAGGAGGCAGCTGCCACTGCTGATCCCA
AAGCAAACAAGAAAAAGATAAGGACAAAAGAGAAAAAGTTTGTGGAACCATGGAATTACTCTACCTCT
AAAAAATGTCAGAAAGAGACGGTTCGGGAAGACAGCAAAGAAGAAGTACATTGAATCTCCGGATGTGGAA
AAAGAAGTCAAGCGTCTGCTGAGTACAGATGCAGAAGCTGTCAGTACTCGTTGGGAGATAATTGCTGAAG
ATGAAACAAGGAGACAGAAAATCAAGGCTGGACATCTCTTCCAGGAATGTCTGGCCACAGGCAAGG
CCATGACTCACTAGAGCATGATGAGCTTCGAGAGATATCAATGACCTCAGCAGCAGCAGTGAAGATGAA
GAAGACGTAACATCCTTGACACTGAGGAAGATCTGGAGAGACAGCTACAAGACAAGCTAAATGAATCCG
ACGAGCAACACCAAGAAAATGAAGGAACCAATCAGTTGGTTATGGGAATTCAGAAGCAGATAGATAACAT
GAAAGGCAAGCTCCAAGAAACCAAGACAGGGCAAAGCGCCAGGAAGATCTCATTATGAAGGTAGAAAAC
CTGGCTCTCAAGAACAGGTTTCAGGCAGTCTGGATGAGCTCAAACAAAAAGAAGACCGGAAAAGGAGC
AGCTTAGCTCTTTCAGGAAGAGCTGGAATCACTTCTGGAGAAGTGAAGTCTGTTGACTTTAGTCCCTTGA
CTAGTCAGCTCAAGAAAAAATTATTTTCAGAGTAAGGTCTTGAATACTTTTTTTTTCTTTGAAATGTTT
TTATCCCTGAGTTTTGTTGAAAATTTTTGATTGTTTTTCAGAACTATGGTGGTGGGTTTTTCACTCTTA
AGAAAATAGGATTACATAGACAGGTTTGTGTAGGGAAGATATAAGAAAGATGAGTGAAATGAAAACACT
TGAACCTTTAATCACAGGACTTTATTACAATTTGTTTCTAAAAGTTTCTAGTCTTTTTTTTTTTTTTT
TGGTGCTAAGATCAGTAATAAGCTAGCCTACAGGCAGTGAAGTTAATTTCTTTTAAAAGAAATATATTTT
TTTTTACCTCAGCAAGTTTCCATGTTACAGAAATGCAAGTTTTTCATAAAGGTTGTTTCTTATACAATAA
AAGCTACTAATTTTCACTGTTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

```



[View online »](#)

<b>Restriction Sites:</b>	RsrII-NotI
<b>ACCN:</b>	NM_011901
<b>Insert Size:</b>	1026 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC029673</a> , <a href="#">AAH29673</a>
<b>RefSeq Size:</b>	1730 bp
<b>RefSeq ORF:</b>	1026 bp
<b>Locus ID:</b>	24074
<b>Cytogenetics:</b>	18 B3
<b>Gene Summary:</b>	Functions as a component of the DNA-binding general transcription factor complex TFIID, a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. Present in both of the previously described TFIID species which either lack or contain TAFII30 (TFIID alpha and TFIID beta respectively). [UniProtKB/Swiss-Prot Function]