

## Product datasheet for **MC204629**

### Naprt (NM\_172607) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Naprt (NM\_172607) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Naprt  
**Synonyms:** 9130210N20Rik; Naprt1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >BC058800

```

CCCCGCGAATGGAGATGGAGCTGGACAGCGAGGGGCGCATGGTGGTGCGGCCGTTGCTCACCGACCTCT
ATCAGGCTACCATGGCGCTGGGCTACTGGCGTGCGGGCCGGGCGTGCGAGGCAGCAGAGTTCGAGCTCTT
CTTCCGCCACTGCCCTTTCGGAGGCTCCTTCGCCCTTTCGCTGGCTTGCAAGACTGTATGCGCTTCTCTT
CGCGCCTTCCGTCTGCGGGATGCAGATGTGCAATTCTGGCTTCGGTACTACCCCGAGACTGACCCTG
CGTTCTTCGAGCACCTTCGAGCCCTTGATTGCTCTGGGGTACCCTGCGGGCTCTGCCAGAGGGTCCCT
CGCCTTCCCTGGTGTGCCGCTTTCGAGGTGTCCGGGCTCTTCTTCTGGTACAGCTGCTAGAGACTCCA
CTGCTGTGTCTGGTCAGCTACGCTAGCCTCGTAGCCACCAATGCTGCTCGGCTACGCTTGATTGCAGGGC
CGGATAAGAACTACTAGAGATGGGTTTTCGGCGAGCTCAGGGTCCCAGCGGGGTTTTCACGGCCTCAAT
TTACAGCTACCTGGGTGGCTTTGACAGCAGCAGCAACACACTTGTGGACAAGTACGAGGTGTACCGGTG
GCAGGGACCCTGGCCACTCCTTTGTCACTTCATTTTCGGGTAGTGAGGTACCCCGGACCCGATGTTGG
CTCCAGCTTCTAGTGAGGGCCCCACAGTGGACCTGCCTGCCAGGGTAAATTTGTGGCTGAAGCGTGTGTG
CCTCTACTTGGGGCTAGAGGAGCAGGAGCCACACCAGGGGAGCGCGCAGCATTTGTGGCCTATGCTCTG
GCCTTTCGCCGGCTTTCAGGGGCTGCTAGATTCTACAGTGTACGCAGGAGTGGCCTTCTCAACTTCC
TGGCAGTAGCCCTGGCACTGGGGGAGCTTGGCTACAGAGCAGTGGGTGTGAGGCTGGACAGCGGTGACCT
GCTTCAGCAGGCTAAAGAGATCCGTGGCATCTTTCGAAGTGTGGTGCCCAATTCAGATGCCCTGGCTA
GAGTCTGTTCCATTGCAGTCAGCAACAACATTGACGAGAGTGAGCTGATGAGGCTAGCCAGAAGGGCA
GTGAGGTGAATGTCATTGGCATTGGCACCAGTGTGGTCACTTGTCCCAAACAGCCTTCTATGGGCTGTGT
CTACAAGCTGGTGTCTGTAGGAGGCCAGCCTCGAATCAAGCTGACAGAGGACCCGAGAAGCAAACACTG
CCAGGGAGCAAGGCTGCCTTCCGGTTCCTAGGCCCGATGGGTCTTTGCTACTGGACCTACTGCAGTTAG
CAGAAGAGCCCCGCCAAGGCTGGGCAGAACTGAGGGTTTGGCCTCGGGGACCCAGGAGCCCTGCAC
TGTGAAGCCAGCCAAAGTGGAGCCATTGCTGCGACTCTACCTGCAACAGGGACAGCTGTGTGAGCCACTG
CCATCTCTGGATGAGTCCAGACGCTTTCGCCAGCAGTCCCTGAGCCTCCTCCGCCCTGCACATAAGCAGC
TGCAGAGCCCTGCAGTGTATCCGGTGCAGCTGTGAGAGAACTTCGGGCACTGGTGGACAGTTTGTAGTGC
TCGAGGTCCCCTGTGAAAAAAGTCATAAATAACACGGCTCTGTTCCCCCAAAAAAAAAAAAAAAAAAAAAA

```

**Restriction Sites:** RsrII-NotI



[View online »](#)

<b>ACCN:</b>	NM_172607
<b>Insert Size:</b>	1617 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">BC058800</a></u> , <u><a href="#">AAH58800</a></u>
<b>RefSeq Size:</b>	1691 bp
<b>RefSeq ORF:</b>	1617 bp
<b>Locus ID:</b>	223646
<b>UniProt ID:</b>	<u><a href="#">Q8CC86</a></u>
<b>Cytogenetics:</b>	15 D3
<b>Gene Summary:</b>	Catalyzes the first step in the biosynthesis of NAD from nicotinic acid, the ATP-dependent synthesis of beta-nicotinate D-ribonucleotide from nicotinate and 5-phospho-D-ribose 1-phosphate. Helps prevent cellular oxidative stress via its role in NAD biosynthesis. [UniProtKB/Swiss-Prot Function]