

## Product datasheet for MC204642

### H2aw (NM\_178218) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** H2aw (NM\_178218) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** H2aw  
**Synonyms:** Hist3; Hist3h2a  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)

**Fully Sequenced ORF:** >BC063781  
 CTTTTGGTTATCAGTGGAAAGATGTCTGGTCGTGGCAAGCAGGGCGGCAAAGCTCGTGCAAAAAGCGAAGTC  
 TCGTTCTTCCCGTGCCGGCCTGCAGTTTCCAGTGGGCCGTGTGCACCGCCTTCTCCGTAAGGGTAATTAT  
 TCGGAGCGGGTTCGGCGCTGGTGCCTCCGGTGTACCTGGCGGCTGTGCTAGAATATCTGACTGCTGAGATCC  
 TGGAAATTGGCTGGTAATGCGGCCCGCGACAATAAGAAGACGCGAATCATCCCGCGCCACCTGCAGTTGGC  
 CATCCGCAACGACGAGGAGCTCAACAAGCTGCTGGGCCGCGTACCATCGCGCAGGGCGGCGTCTGCC  
 AACATCCAGGCCGTGCTACTGCCAAGAAGACCGAGAGCCACCACAAGGCCAAGGGCAAGTGAGGCCCT  
 ACCGGTCGTTAGTGCAGCTTTTACAGACACCAAAGGCTCTTTTACAGGCCACCCACAACCTCATTAGAAA  
 GCGCTGTACTGAGCCTGTGCGCAGTTACTGGCCAGGCCAGTTCTCCTGACTCCTTCCGTGCATTCTG  
 CATCCCTGTGCGGCGCTTGGGGAGCTCAAGGGAGACCAACGGCTTCATCTGCGGGTGGCCTGCGGGCCGT  
 TGGAGGCGTTCGAGTTTCTCCCCACCCTCCGTGTGAGACACAAGTTGCAGTTAACGTTCTGTGTTTGA  
 CATTAAGGGGATACTTAATAATAGAAATTAGTATAAGGTTCTGGTTTAGATGCTAGTGTCTATTGAAA  
 CTGGGTTAAATTGACATGTTCAAGTCTTCCAGTTAGGACGGGTTTTCAAAGGTTTCTGAGGGGAGACTT  
 GGTATTTCAAATCTACTGTTTATATAATGTTCTGTTACAAGGCGATTTACAACCTTTGGAGTTTGGCTG  
 GTTTCCTTTGGGGGTGGGGAGCCGCTCTAGGCTGCTTACAGCTGCTTAGTGGTAGGATTACAGATG  
 AGCATGGTGGTCTGCAGGCCTTTGGAAGTGGTAAATCGGATCCACTTAAGTAGCTTAACCTTGGGGAGC  
 AAAGTGATTTGTCAGGTTCTGGTAAACTAGGTGTGAGAGTACAGCAGTGCCTTTGTATAAAGTTTAA  
 AGTACAGGTAATCAGGAGCGTGGTGGGAAAGGGAGAGGGGGTTAAATGAGTTACCTTTAAATCG  
 TTTTAAATGGCACAGAATTGGGCCTCTCCCTGTTAACTGCTGTACTACTTTAAAAATTTAACAATACAG  
 GCCACGAAGACAGTTATATAATGTTTATTATATAATTTATATAGACTGTACTGAAACTGGATCTGTTCAA  
 GCCTGGGCTCCTCAGACCCGTCTGTAACACAGAGCTGCTATAAAATACCAACGTCTTAAAGCCATTT  
 GTCTAAGGATAAGCAATTTCTTAAACGGTGAGTCTTAAACAGAGTGATCAATAGACTGCGTATAGATTAT  
 ACATCTGCAAATGGGTGGTCTCAGTAAAACAAAGCCAGAGGGTTGGACCTCCAAAAAAAAAAAAAAAA

**Restriction Sites:** RsrII-NotI



[View online »](#)

<b>ACCN:</b>	NM_178218
<b>Insert Size:</b>	393 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC063781</a> , <a href="#">AAH63781</a>
<b>RefSeq Size:</b>	1607 bp
<b>RefSeq ORF:</b>	393 bp
<b>Locus ID:</b>	319162
<b>UniProt ID:</b>	<a href="#">Q8BFU2</a>
<b>Cytogenetics:</b>	11 B1.3
<b>Gene Summary:</b>	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H2A family. [provided by RefSeq, Aug 2015]