

Product datasheet for MC204759

Nufip1 (NM_013745) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nufip1 (NM_013745) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Nufip1
Synonyms:	Nufip
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC056192

```
GACAAGAGGTAATCCCTTGTAGGTGGCATGGCTGAACCGGCACCCGCCGTGTGGCCTTCGGCCCCGGAT
CTGACTCCCGCGCCCGGACTCCGAGCGAAGCTGCGCCGCGCGGACAACCTGGGTGTACTGGGCCATGC
TGCCGCCGCCACCACTCCTCTGTCTCCAGTGGCGGGCTCGGAGCAGAGCCGGAAGGACAGCCTCA
CGTGCTGCCCCAGCCTCCTCCGGAGCGCTTCCACCGTTCGATGCTCAGATTCTCCCGCGCGCAGCCT
CCTTTTCGATGCCAGGCTCCGCCGATGCTCAATCTCAGTTCAGCGGCCAGCAGGCTGGAATTTGCAAG
CCTCCACGCCCTTGGTACTGGGATTGTCTCCTAATGGTTTTTCCACGTACCACACGTCTTACCAATCTCC
GGTTACACATTCTATTTTCCACGATCACATGATGCAAAATTCATTTGCCTCAAAAACAGAAAACAGAAA
ACGAAGAAAAGAAAGGAACCGGTTTTCCATTTTTCTGTGATACCTGTGATCGTGGCTTTAAAAATCAAG
AAAAGTATGACACACACATGTCTGAACATACAAAATGCCCTGAAGTAGATTGCTCCTTTAGTGCACATGA
GAAGATTGTCCAGTTCATTGGAGAAATATGCATGCTCCTGGTATGAAGAAGATCAAGCTAGACACTCCA
GAGGAGATTGCAAGATGGAGGGAAGAGAGACGAAAAAACTACCCAACCTCTGGCCAATATTGAAAGGAAGA
AAAAGTTACAACCTGGAAGGCAAGAGAGAGAGAGTGTGACAACCACACAGTATGGCAAAATGAAGGG
CATGTCCAGACATTACAGATGGCAAAGATCAGAAGTCCCTGGTAAACATCATAAATGGAGAAGAGGTGGG
GCTAGGCAGAGAGCAGTCGTTGGACTAGGAAATCATGCCAGAGATTGCAAGCCTGAAGTTCCATCCAAGG
CTAATGTGGACCTCTAGGTGCTTTGATCCACAGCGACTCTGAGTCTGATAAGGAAGAGAAAAGCGCAAAAG
GACTGTAGTGCCCAAGGAAGTGACACCAGCCTTGTTGCTGATGAGCAGCTACGGCGATGTCTCGGGG
TCAGACAGTGAGCCAGAAGAAGTCCCATCAAGACTGAAGCAGAAGTTCTGGCAGAAAACCATGTTCTTC
ATAGCAGTCCCTCCAAGAGTCCAAAACAAAATGTTCAAACAACCTGGTAGAACAGTCTCAAGATCCAAGT
GGAGAACCAGAGGAACGACTCAGAAAAGATAAGCCTTAAGCGGAAAAAAGGCCACTGTCAACCCATTATTT
GAACCAAGAACACATCACCCATATCTCCTGGAATGCTTCTAGCTCCAGATATTGACATGAAAGAAATG
TGATTTTGCAGTGTGTTGATATATCATCAAAAAAGACTTTTTTGGACTTAATACTGATTCTGTGAAAAAC
TGAAGAAGTGTAGATGTCTCCGTTTTCCACACTCATCTGTAAGCACCTGAGAGAACCACATCCTTGAGT
AAGGAGAGAAAATCAGCAACTTCTGCAGCTGGGATCAGTAACTGAACTGGAAGCCTCATGGGCGTCATGT
TGGATTTAAGTCTGCTCCTCAGGCTTTATGCAATAAATGCATAGCTGGTCTTTTTTAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_013745
Insert Size:	1455 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC056192 , AAH56192
RefSeq Size:	1719 bp
RefSeq ORF:	1455 bp
Locus ID:	27275
UniProt ID:	Q9QXX8
Cytogenetics:	14 D3
Gene Summary:	Binds RNA.[UniProtKB/Swiss-Prot Function]