

Product datasheet for **MC204792**

AI467606 (NM_178901) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AI467606 (NM_178901) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	AI467606
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC064101

```

CCAAGATTGAAGCAAGGTGACAGGAAACGGGATATCCTCTGGTTGCTGTCCTTTCTCCAGTGGAG
ACGTTCCCTGACTGTTCTGAGTCCACACTGAATCTCCAGATGCCAGTGACGCCTCAGCAGCCCTGGG
CACACGGAGGGGCTCCCTGAGCCTACAGCGGAGGCAGCAGTGTGGGTAGTAATCCCTGTGGACCTGTGTA
TTCCCATCATGCTGGCCTGGCCTCTCTCACAGCCTTCTCATATAACTACTGCTGTAAGTGGCCGAACG
CCTGTTCCGCCGCCGCAACAGACCCAGCCAGCGTGCCCCACCCTAGTCTGGCGCCCTGGAGGAGAA
CTGTGGATTGAGCCCACGAGCAGCGCCCGGAGCTCTCTGAGGATTGGTACGGCTCCTCCATGCCTCTGC
TGATGGACAGGGCTCCAGGTCCCCAACACCCGGTGGCACCTTAGAGGGCAGAGCAACTGCTCCACCAGC
CACTTCAGCCCCATACTTTCCCTCAGTTCTCTAGTTCCTCAGACCCACCAGAAGTCCCAGCCCAGAGC
ACCTTCTGGAGGCCTCAGACCCAAGAGGAGAGGCCCATGACACAAGCCTGGTGAAGTGGGTTGGGCTG
AGCCCATGCCAGAGGCTGGCTTGCAGTTGGGAGCCCCAGGCCCTGGAGGCCGAGGCAAGGGAGCCGGA
ACCTGATTGGGGTCTCCAGCCTCGGGTCAACCCTGGAGCAGATCTCAGCATTCTGGAAGAGGGAAGGCCGG
ACCAGCGTGGGCTTCTGAGTCTTCAGGACCTGTAAGACTAACTCCACCTCAGAGGGACTACAGGGGTCA
GCTCCAGCCTTAGGCCTGACCCTTGCCTGAAATTTTGTGTGAAATTTCTGTCTTCTGCACCTGGCTG
GTGGGCAGGCTAAGAGTGTGGACAGCCTGCACTCCAAGGCTTGCCTGAAGTTCAGCTTAGGTTAAATGC
CTAAGTCGTGACTCTGAGAATTGAAGGCAAGTCTAGCTGGCTACAGGCCAGAAGGCTTCAAGATGGCAG
TCAGGGTACAAACTTACCTAAGACAGCCTATCCTGTATGGGATGCTGACCATCAGCCAGAGCTGGTAGT
TATGAGTGTGAGTGTGGTCTGGAGATGGCAAGAATCTTGCTCTCAGAGCCACAAACCTATCTGTG
CCTGGCCACTCTCTCCACCAGTGTGAGGCCCTCCAGGAAAACACAGAGGACAAATCCAGGAGGTGTCTT
TGATTCCGTACAATAATGGCCTTGGCCACTAAACAGCAACATGGGATTTGATGGCCCACTGGAGAGG
CTGTCCTTGCTGGAATTTGGACACAGAACCAAAACAGATCAGTAGGCTCCTCTGCCCTGACTTATGG
CTGCACACATCCTACAGTACAGTTTTATTCTGTGTTAAGACAAGCATGGCCTATCGAGAAGTGAAGT
TGTTCCAATATGGGAATGAAGATTTTTTCCCAGAGAAAGGTCCTCCCAAAACATAGCCTCTGCCATCT
TTCCCTCCAGCTTGTGGAGTGGCTAGAGAGGCAGGTTTGTGTGGGCAGACTTACTTACAGTATCC
ATACCATAGAAACGGAGTCCCAGTACCCGTGTGCCACCAGCTGACAACCTCTTAGGAGAGGGGTTCC
CAGTGGGAAGATGGAGCTTTCCGCATCTGAAATCTGACTTTCCAAACTCAGTCAAGAGTATGACGAT
CCTGCCGTATGTTACTTTTTGTTTCATTATTTACTTTATATTTTACTTAAATCGGTTCACTGTATGCT
TAATAAATTTGTTTTAAAGTGAAAAAAAAAAAAAAAAAAAAAAAAA

```



[View online »](#)

Restriction Sites:	EcoRI-NotI
ACCN:	NM_178901
Insert Size:	678 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC064101</u> , <u>AAH64101</u>
RefSeq Size:	1863 bp
RefSeq ORF:	678 bp
Locus ID:	101602
UniProt ID:	<u>Q8C708</u>
Cytogenetics:	7 F3