

Product datasheet for **MC204840**

Crct1 (NM_028798) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Crct1 (NM_028798) Mouse Untagged Clone
Tag: Tag Free
Symbol: Crct1
Synonyms: 2300002G24Rik; AA589624; Nice-1
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC060280
 CAAGTTCCAGTTGGAGAAGTCTGAGTCTTCTGTTGGAGCCAAGGTCTAATTCATCTTGACAAACGCCG
 AGATGTCTCAACAGGGCCCAAAGGCTTTTCCAAAGGGTCTTCCAGAGTTCCGCCCCCTGCCCGCCCC
 GACTCCTGCCCGCGTCTCCTCCTCCTTGTCTGTGGCGGAGGCTGCTGCAGGGGCTGCGGGGGCGGC
 TGCGGTGGTACTCTGGCTGCTGCGGCTCCAGTCCACCGGCTGCTGCTGCTTTCCAGGAGACGCCGCC
 GACAGCGCTCTGGCGGCTGCGGCTGCTGCTGCGGAGGAAGCCAAAGGTCTCAACGCTCCTGCAACAATCA
 GAGTTCAGGCTGCTGCTCCGGTGGCTGCTGAAAGCCTACAGAGCTTCTCCGGAACAGCTGGCATCGACA
 GCGTCTTTGCTGCTCTACCCAGTCTCCATAGTAGTAACAGGAAGGACCGGCCTAGGAGTTTGGTTTCC
 CTTAAGCTAACTGCACTTTGATGTTTCAGAACTTACCCTGTTCTTAGACCTGCTCCTAACGTTGATTT
 TCCTCCAGGGGATACAAGTCAGGGGTGGGACTTCAGATCAGCTATGGGCAGGAGTTGTGCACAGCAG
 TGTTCAAAACAGGCCTCCTGATTTGATTGTCGCCTGGAGACATAATAATAAACTTCTACCTCCTAAAAA
 AA

Restriction Sites: RsrII-NotI

ACCN: NM_028798

Insert Size: 309 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC060280](#), [AAH60280](#)

RefSeq Size: 751 bp

RefSeq ORF: 309 bp

Locus ID: 74175

Cytogenetics: 3 F1