

## Product datasheet for **MC204899**

### **Dbil5 (NM\_021294) Mouse Untagged Clone**

#### **Product data:**

**Product Type:** Expression Plasmids  
**Product Name:** Dbil5 (NM\_021294) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Dbil5  
**Synonyms:** ELP  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)

**Fully Sequenced ORF:** >BC048474  
AGCCGCTGCCGTTCAAGGCGCACGCGTACTCCTCTCCAACAAGCTCGCAGCGTCGAAGCCACCCTAGCCA  
CAGCGCCCCCTCTCGGCCTAGGCTTCCCAGCCTTCTTAACCACAGTACCCCGCATCACTTAAAAACAGAA  
GCATTGCCTCTGGTATACCAATGAGCCAAGTGGAGTTTGAAATGGCCTGTGCTTCACTCAAGCAGCTGAA  
GGGCCTGTGAGTGATCAGGAGAACTGCTGGTGTACAGCTTTTACAAACAGGCCACCCAGGGCGACTGT  
AACATCCCTGTCCCTCCTGCCACAGATGTGAGAGCAAAGGCCAAATATGAGGCCTGGATGGTAAACAAAG  
GGATGTCCAAGATGGATGCCATGAGGATCTACATTGCTAAAGTGGAAAGAGCTGAAGAAAAAGGAGCCATG  
CTAAGAGCATTTTGGCAGACACAAGGAAGTAGCATGGCTCCACTGGTAGGGGACGGGTTTAAAGATTCCTG  
ATGGCTGAAACACACTGAAAGGGTAGCAGCGTTAGCAGAAACACCAATAAATCATTCAAACACTACAAAAA  
AAAAAAAAAAAAAAAAAAAAAAAAA

**Restriction Sites:** Sfil-Sfil

**ACCN:** NM\_021294

**Insert Size:** 264 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC048474](#), [AAH48474](#)

**RefSeq Size:** 584 bp

**RefSeq ORF:** 264 bp

**Locus ID:** 13168

**UniProt ID:** [O09035](#)

**Cytogenetics:** 11 B5

**Gene Summary:** May be involved in the energy metabolism of the mature sperm.[UniProtKB/Swiss-Prot Function]