

Product datasheet for **MC204902**

1700016D06Rik (NM_024271) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700016D06Rik (NM_024271) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1700016D06Rik
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC048664

```
GTATCCGAGAGATCAAGGCTGGGAAGCTGAAGACCTCCAGGGAAGGCCTCCTCCTTCTGCTGGAGGG  
ACCATGACAGCGAATAGAGGAGCAATGACCTTCTTATCTGTCTTCTGCTGTGCTTTGGCAAGACGTAC  
AGCCCTGGCCGAATAGAGAGCTTTCAGGTTCAAGTGGTAAGCAATGGCCAAATTTAGAACTCTGGATGA  
TGACATAGCTACAGTGTTTGATGAAATTCTAGTCCGGGATATTCTGGAGCCAGGGAAGCACCATATTTT  
GAAATCAAAGCCCATCAACAATTCACAACAAAAACAACGAAGGAAAAAGAAGTCCACACAAAAGAAC  
CCATACCTTCTAAGACATATCAGAAGCAGTCATCTTCTGGGACCATGCACACCCCTTCTTTGACGATAA  
GCAGAAAGAAATACTTTTTCAGCTAAGAAGCCTTGAAGCTCTGGAGAAAAATGATTGACAAGATTCGAAGA  
GCTATTGAAACTAAGCTCAAGAGAAGACAGAAGCTCCAGCAAAGTAGAAGCAGTCAACTCCTGGGGAAGC  
CAATCCAGCCCTAAGTACTAGGCCTGGCGAGAGGGACAAGCCACTGCCGCACATGGGGACAGAAGAGAGT  
GGAAGGAGCCATCGACCTGTGTGAAGTACTCAATGAGCACCAATAAAATTAGCATATTTAAATAGTA  
AAAAAAAAAAAAAAAAAAAAAAAAAAAAA
```

Restriction Sites: Sfil-Sfil
ACCN: NM_024271
Insert Size: 501 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC048664](#), [AAH48664](#)

RefSeq Size: 728 bp

RefSeq ORF: 501 bp

Locus ID: 76413

Cytogenetics: 8 A1.1