

Product datasheet for **MC204912**

Hrasls (NM_013751) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hrasls (NM_013751) Mouse Untagged Clone
Tag: Tag Free
Symbol: Hrasls
Synonyms: 2810012B06Rik; A-C1; Hrasrs; HRSL1
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC048482
 AGGGCGATGCCGATGGAGTCACCCAGGGTCTGCTAGCGGCACCGCAGCGAGCCCGTGGAGAACCCGGGAA
 GCCCGTCACTGCCTCCAGGGCACTAAGACAACCCAGCTTGAGCAGGGAGATGGCGGTAATGATTGCTTC
 AGTCTGACCTATCCTCACAAACCCACACCCAGGAGACTTGATTGAAGTGTTCCGTCCTTGCTATCAGCACT
 GGGCACTGTACTTGGGTGATGGCTACGTGATCAACATTGCACCTATAGATGGCATTTCGCTCATCATTAC
 AAGTGCTAAGTCCGTGTTGAGCACAAGGCCTTGGTAAAATGCAGCTTTTGAAGGATGTTGTGGGAAAT
 GACACATACAGAATAAATAACAAGTACGACACAACATACCCTCCTCTTCCGTGGAGGAGGTGATACAAC
 GGTGAGAGTTCGCTATTGGGCAGGAAGTAGCCTATGACTTGCTGGTCAACAACACTGTGAGCATTGTTGTAAC
 CTTGCTGCGCTATGGAGAAGGAGTGTGAGAGCAGGCCAACCGAGCAATCGGCACCATCGGATTGGTGGCA
 GCTGGTATTGATATCTTCACATTCCTCGGCTTGTTCCTCCAAAAGACAAAAGAAATATTAGCAGTTTA
 TTGAAGAGATTGGAATGAAGAAATTTGTGAGGAGAAAAAAAATCCTAGGGTGAATACTTATTTTGAATG
 CATCATTATTGCTCATGGTCCCATGATGGATGGCAGACTCTGTAATAAATTGCTTGTGATAAAAAAAA
 AAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: SfiI-SfiI

ACCN: NM_013751

Insert Size: 504 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC048482](#), [AAH48482](#)

RefSeq Size: 792 bp

RefSeq ORF: 504 bp

Locus ID: 27281

UniProt ID: [Q9QZU4](#)

Cytogenetics: 16 B2

Gene Summary: Exhibits both phospholipase A1/2 and acyltransferase activities (By similarity). Shows phospholipase A1 (PLA1) and A2 (PLA2) activity, catalyzing the calcium-independent release of fatty acids from the sn-1 or sn-2 position of glycerophospholipids (PubMed:21880860). Shows O-acyltransferase activity, catalyzing the transfer of a fatty acyl group from glycerophospholipid to the hydroxyl group of lysophospholipid (By similarity). Shows N-acyltransferase activity, catalyzing the calcium-independent transfer of a fatty acyl group at the sn-1 position of phosphatidylcholine (PC) and other glycerophospholipids to the primary amine of phosphatidylethanolamine (PE), forming N-acylphosphatidylethanolamine (NAPE), which serves as precursor for N-acylethanolamines (NAEs) (By similarity).[UniProtKB/Swiss-Prot Function]