

Product datasheet for **MC204923**

1700129C05Rik (NM_026461) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700129C05Rik (NM_026461) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1700129C05Rik
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC049552

```
AGTCAGTAACTGCCATCAGTAAAGTAAGAGCTCTGCAGGGTCTGGGATGATGAGAAGGAAGCGGCCAAGG
ATTTTGGGAAGTCTACGCAATCCTGTTCCCCCTAAGGATCTTTCTGCCCCTCTGATCAACCAAAAAACA
TTCATTTCTTCCACCTCAGACTCCCATGCCTCAGCCTACCTTTTCGTACCAACCCCTTTTCGTTTTCTA
CCCTACCTCCTTCCCTGCCCTGAGGAACACTTCTACTTCCCCACCTTAATGATTGCTATGCCAGCTCC
TCCTGTGTAACAAACTTTCCTCTAAGCACCCCTCTTATTTCCAGACCCCTTCTTCTGTCTACCCCTCA
TGTCTACCTCTTCTCCAAGAGCACTGCTGGATAGTTGCCTAGTTCCTCTTTCCCATCTGGAACTAGCTA
CTTTCAGCCAGAGGAGCAAATCCTTCAGGAACTGAACAAGATTTCTCTAGACCCAGAGCTTAAAATTGGA
GGAAGCAGATATGACCTCACATCCTTGGTACCAGGGAATCTGCAGCCAGGGACAGAGACTGGAAATCATG
AAGTGGAAAGATGCTCAATCCCTGGATGATCTTGATGCCTTCACTCAGGAAACACTGGACAGAATGAATTG
TCATTTGGTCAAGAACGCAAGTTTACCAAATAACATTTTCGGGATGGAGAAACCAGAGGGGAAACTGAAG
TGACAGGACCTCAGGATTTCTAAATAACTTCAGGACCGAAAATAATGCCACTATCCAAAAAGCTAAAATA
TGATATCTGGGAAAGATGGTCATCAGAAGCAGCATTATCCCATGTAGCCATTCCATGGATAAAATAAAC
AAAGTGGACATGGTTCACAGAAATGTTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AGAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
```

Restriction Sites: Sfil-Sfil
ACCN: NM_026461
Insert Size: 657 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC049552](#), [AAH49552](#)

RefSeq Size: 942 bp

RefSeq ORF: 657 bp

Locus ID: 67932

Cytogenetics: 14 C3