

Product datasheet for **MC204926**

Abhd16b (NM_183181) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Abhd16b (NM_183181) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Abhd16b
Synonyms:	1700095H12
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC050777

```

GCCTCCGCAGACACCGGCTGGGCCAGGCCCTGCATCTCCAGCTGTCATTTACAACACTAGGACTGGGCCTTC
AGGGAACAGCTCAGATCTTCAGTGCCCATGTGTGCATCTGCTTTGTGAAGGCAGTGTGCACGTGTTCA
AGATCTACTTGACAGCAAACACTACAGTTACAACCTCCGCAGCTGGCCGGTGGACTTCCGCTGGGATGACCT
GCATGCCCAAGCACTGGTAATAGCAGTCAACGTGCGTTGACCTGTGCTGCAGCTGCTGCAGGTGTGTGG
CTGCTACACGATGCTGCACTGGGAGGCGACACACTCACAAGGCCTCCACGTGGAGCTCGGAGCCAAGTGC
AGTGCCTCCTCCAGCAGATCCGTGAGCTGCCAGCCAGCTTGCTAGCTATGCACTGGCCCACTACTGGG
CCGCTGGCTTGATACCCCGTTCTATGTTCTTGATGACGCGTGCATTAATGCCTCTACTACAGCAGGGC
CAGGAGCGACTAGTGGACCGCTACCGTGGGCGCGTGCCAAGCTGGTGGCCTGCGACGGCAATGAGATTG
ACACCATGTTTCATGGATCGCCGTCAGCATCCCGGCAGCCATGGCCGCGGCTTTGTCTGGTCATATGCTG
TGAAGGAATGCTGGCTTCTATGAGATGGGCTGCCTGTCCGCACCTCTGGAGGCTGGATACTCAGTGCTA
GGCTGGAACCACCTGGCTTTGGAGGCAGCACAGGTGCACCTTTCCACACAACATGATGCAAATGCAATGG
ATGTAGTGGTCAAATATGCGCTGCACCGTTTGAACCTCCCGCCAGCACACGTGGTGGTTTATGGCTGGTC
TATAGGAGGCTTACAGCCACTTGGGCCACCATGACTTACCCAGAAGTGGGTGCACTAGTGTGGACGCC
ACCTTTGATGATCTGGTGCCACTAGCGCTGAAGGTTATGCCCCAGAGTTGGAAAGGGCTGGTAGTACGCA
CTGTACGGGAACACTTCAATCTCAATGTTGCTGAGCAGCTGTGCTGCTACCCAGGCCCACTCTTGTCT
CAGACGTACACAAGATGATGTTGTCAGCACCTCGAACCATATACCCTCGATGCCATCTTGCCAAGTGGAG
GGGGATGTGGAAGGTAACCGTGGCAATGAGCTGCTGTTGCGCCTTCTACAGCACAGATATCCCTCTGTGA
TGGCTCGCGAGGGCCGCACCGTGGTAACCCGCTGGCTACGTGCCAGCAACTTGGCTCAAGAAACTGCATT
ATATGCACGCTACCGCGTTGATGATGAGTGGTGTAGCAACACTTCGTTCTTATCGTGAGAGATGCCAG
AAGGAGCTGGACGATGCCGAGGCCCTGGGGCCACACGGACTTTCTTTCTTGGTTCGTGGGCCAAGGTC
TGAGTGCCCGGCGCGCCGCGCAGCTTGCCTTGGCACGCAGGCACCTTAAGAACCTGGAGGCAAC
TCACTGCAGTCCACTGGAGCCTGAGGACTTCCAATTGCCCTGGAGACTATAGCCACTGTTACTGTGACTC
ATGCATTTGGGCATCTGTACACGTTCTAGCTCCCCCAACTGTTCTCCTTACATCTATACAGAAAAG
CTGGCCGTGCTGGGGATCTTGGTGTCTGGGAGAGTAGGCAGTATATCTTCTTCACTTTCTTCTGTTCT
TTCTGCTTCCATCCTTTCTTCTTCTTCCAAACATCTGAATGTCTGTGCCTAATTCTTAAAAGCAATGTA
AAGATGACCTGATCCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

```



[View online »](#)

Restriction Sites:	Sfil-Sfil
ACCN:	NM_183181
Insert Size:	1425 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC050777</u> , <u>AAH50777</u>
RefSeq Size:	1796 bp
RefSeq ORF:	1425 bp
Locus ID:	241850
UniProt ID:	<u>Q80YU0</u>
Cytogenetics:	2 H4