

Product datasheet for **MC204968**

1700001C02Rik (NM_029285) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700001C02Rik (NM_029285) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1700001C02Rik
Synonyms: 1700047H18Rik
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC060996
 GAGCCGAGCCACCGCGCCGGGGCTCGCAGACCCACCATGGCCTTCCGTAGCGCGGGCACACTAATGACC
 GAATTTAATGCCGCTTCGTGCCCCAGCTCTCATGCCCGGTACAAAGGCCATGTTCCCGGTGTGGCTT
 CCTCCTTTGGCTCTTCTTATGGCAACACCACCTTCAAGTACTTCCAGGACCTCCGCAACACAGGCTTGGA
 GAAGAGCTATGCGCTCCTCAGTGGGGTTGCTTTCCACCATCTTCTCCCCTAACTCCCCCTGGTTCTG
 ACTGATAACTCACAGAACTGGGATCGATGGCTGCACCAGCCACCTACAGCCGATCCAACCAGGACAACG
 GGCGCACCACTGAGCTCACTAACTTTTACAAGACAGTCCAGCAGCAGCGGAAGTATTACCAAGACAAGAC
 TGGCATTCTGTCCCGTGTTCCTTACTTTGTACTTCTGTGAAGGAGCCGGACCGGTACCCCTACCCACC
 GACCTGCCTCCTCTGAGTGAGAAGAGTAAGTGAACCTACTAAGATTATCCCTGACTACAAGAGGACCT
 ACCAGACCTTCCCATCAGGGAAGAGAGTTTCTCTCAAGAGCGCCAAAGCGGGACCTTTACTTCGAGTT
 CAGAGCCTGAAACCAGGTGGTTGACCCAGGGAACCTTGGCCTGACCAGGCTGGCTTGCTTGAATAAAA
 TCTTCAGCCACAAGCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: Sfil-Sfil

ACCN: NM_029285

Insert Size: 603 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC060996](#), [AAH60996](#)

RefSeq Size: 745 bp

RefSeq ORF: 603 bp

Locus ID: 75434

UniProt ID: [Q9DAS2](#)

Cytogenetics: 5 B1