

Product datasheet for MC204975

Cyct (NM_009989) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cyct (NM_009989) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cyct
Synonyms:	T-Cc
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC059728 AGGATGGAGCCCCGCCATCAGCTGGTGCTGGGTCCTGCCTTTAAGAGGAGTCCTTGGGCACAGCAGTTGCA CCTTACATCTTGGCCTCCCTGCATCAGTGGACCAAGGTGGCCAGGGCACGGCTGCTGTGATTGTGAAT ACAAAGTATAAGGAAACAAAGTCCAGTTTAAGGTCTTATTGCAGAGCAAAATGGGAGATGCTGAAGCAGG CAAGAAGATCTTTGTTTCAGAAATGTGCTCAGTGCCACACGGTGGAAAAAGCGGGAAACACAAGACAGGT CCAAATCTCTGGGGCCTCTTTGGTCGAAAGACAGGACAAGCACCAGGATTTTCTTACACGGATGCAACA AGAACAAAGGTGTTATCTGGTCAGAAGAACTTTGATGGAATATTTGGAGAACCCGAAGAAATATATCCC TGGCACTAAAACGATCTTTGCTGGTATCAAAAAGAAGAGTGAGAGAGAAGATCTCATTAAAGTATTTAAA CAGGCCACATCCTCATGAAGTGAAGTAGTATTTTCAGCTTCAAAATATTCGCTACCAAGTTTCAAATATTGT GTTTAATAAGAAGCATATCTGATCTTTTTGATCAGCCTGAAGTAATTTAACTTGTTAAATAAAGGTTG ATATTTTGTATTGTTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
Restriction Sites:	Sfil-Sfil
ACCN:	NM_009989
Insert Size:	318 bp


[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC059728](#), [AAH59728](#)

RefSeq Size: 676 bp

RefSeq ORF: 318 bp

Locus ID: 13067

UniProt ID: [P00015](#)

Cytogenetics: 2 44.94 cM

Gene Summary: Electron carrier protein. The oxidized form of the cytochrome c heme group can accept an electron from the heme group of the cytochrome c1 subunit of cytochrome reductase. Cytochrome c then transfers this electron to the cytochrome oxidase complex, the final protein carrier in the mitochondrial electron-transport chain.[UniProtKB/Swiss-Prot Function]