

Product datasheet for **MC205000**

Stard6 (NM_029019) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Stard6
Synonyms:	170001IK09Rik; 4833424I06Rik; 4933429L05Rik; AI644424
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC061022 GGCACTAAACTCTGGACTAAATTAGATAAAGGTAGATGGACTATAAGGCAATTGCCCAACAACTGCTG AACAAAGTTTTAGCTTACAATCAAGACCTATCTGGCTGGAACTCATTAAATCTTCAAAAAAGGTAACGT TTCCAGCAAGACCTCTAGAATATTCACGGAAACCTATACCGCTTGAAGGAATAATCCCAGAATCAGCA GCTCATCTCTGACTTCCTCTTCAAACATGATCACAGAGTCTCATGGGACAAATCATTGAAAGGGTTCA ATGTGATACACAAGATTGATTCGGACACACTCATATGCCACACCATTACACAAAGTTTTGCCATGGGCTC AATTTCCCTAGAGACTTTATCGACTTAGTACACATCAAGCACTATGAAAGGAACGTGGACATTATCAGT ACTAAAAGTGTGGACTTCCCAGGATATGCTCCAACCTTCGACTTATATCCGTGGCTTTAATCATCCCTCTG GCTATGTGTGTTCCCTCTGAAGGAAAACCCAGCCTATTCCAAGCTAGTGATTTTTGTCCAGACAGAGAT GAAAGGAAAATTGCCTGCATCGGTAATTGAAAAGTCTATGCCTTCCAACCTAGTAAGCTTCCTCTCAAT GTGAAAGATGGCGTAAAGACGTACAGAATCCACCCATACGTGCCCGGCACAGTTCGCACTCCTCTGTGC ACAAAAAGAAAGAAGGACATTCTGCCATTAACCGTAGCGATTTCAGACAGCACAGCCTACTAGTAAGTG GAGTTTGTCTGCAGAGTACTCCTAAGTCAACGTGAATAAGCCTGTTCAAGTGTCTCTAGGCATTAACCT TGTACAGAGTACGTAATAATTTGTCAAAGTAAAGACTTCACTGGACAAGAAAAGTCAATTAATACAG TTTTTAAGCCAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
Restriction Sites:	Sfil-Sfil
ACCN:	NM_029019
Insert Size:	702 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC061022 , AAH61022
RefSeq Size:	950 bp
RefSeq ORF:	702 bp
Locus ID:	170461
UniProt ID:	P59096
Cytogenetics:	18 44.44 cM
Gene Summary:	May be involved in the intracellular transport of sterols or other lipids. May bind cholesterol or other sterols (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript. Variants 1, 2, 3, and 4 encode the same protein.