

## Product datasheet for **MC205010**

### Nkain4 (NM\_021426) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Nkain4 (NM\_021426) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Nkain4  
**Synonyms:** AB030182; C030019F02Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >BC061107  
 ATACGCGGAGGATTCCAGCCCACTCTGGCAGCAAAGGCTCCCACGGGCGAAGTGAATACCATGGGCTTTT  
 GCTCCGGCCGCTGCACGCTTCTTGGCCTCTGCGCTCTGCAGCTGGTCACAGCTCTGGAACGGCAAGTCTT  
 TGACTTCCTGGGTACCAGTGGGCACCCATCCTGGCCAACCTCACCCACATCATCGTGGTATCTTGGGG  
 CTCTTCGGCACCAATCAATACCGGCCACGCTACATCGTGGTGTATGTTGTATGGGCGAGCTGTCTGGGTTA  
 CCTGGAATGTCTTCATTATCTGCTTCTACCTGGAAGTGGGGGGCCTCTCAAAGGACAGTGAACCTTTTAC  
 CTTCAATCTCTCTGGGCATCGTTCCTGGTGGGAAGAGCATGGGCCAGGCTGTCTTACGAGGAGGAGGCC  
 ACAGCTGGCCTGGGTGCACTGCATGGACAGTCTTGGTGGTGGTGTGGCTGTGCCATGGTGCACAGCT  
 ATGTGGAAGCCCTGCATAGTGGCCTGCAGATCCTGCTAGCGCTCCTGGGTTTTGTCTATGGTTGCTACGT  
 GGTGAGCGTGTGCACAGAAGAAGAGGACAGCTTTGATTTTCATTGGTGGATTTGACCCATCCCCCTGTAC  
 CATGTCAATGAAAAGCCGTCCAGCCTCTTGTCCAAGCAGGCATATTTGCCTGCATAAGTGAGGGGATGGC  
 TGACCCTGCTCCTGCAACCCAGCCTCAGCCACGGAACAATGAGGAGACCCAGTCTGGAGTAGCAGCT  
 CTGCTTCCTGCCTGCCTAACATTTTTTCCAAGTTCACCTTGGATACCAGTAAGGGGCGGGGACTGC  
 AGGACAGGGTGGCTACAACCTGTGAGTTTGAATATCTTTTCTTTTTTTTTTCTAAAAAAAATACTATG  
 AAAATGCAAGAGAGACTATGGAACTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

**Restriction Sites:** Sfil-Sfil

**ACCN:** NM\_021426

**Insert Size:** 627 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC061107</a></u> , <u><a href="#">AAH61107</a></u>
<b>RefSeq Size:</b>	966 bp
<b>RefSeq ORF:</b>	627 bp
<b>Locus ID:</b>	58237
<b>UniProt ID:</b>	<u><a href="#">Q9JMG4</a></u>
<b>Cytogenetics:</b>	2 H4