

Product datasheet for **MC205023**

Gtf3c6 (NM_026113) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gtf3c6 (NM_026113) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gtf3c6
Synonyms:	2410016F19Rik; AU019813
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC048446 GCGTGGCCAGTGATGGCGTCTCCGAACGCTATGGCGGCGGCAGTGGATCCCATTGGGAGGACTTGGAGG AGGAGGAAGAGGAGCTAGAGGAGGAGCAGTTTGTCTGGTAGAATTGTCAGGAATCATTGATTCAGACTT TCTCTCCAAATGTGAAAATAAATGCAAGATTTTGGGAATTGACACGGAGAGGCCCATCATGCAAGTGGAC AGCTATGTGTTTCTGGGAGTATGAAGATACGCTTGGTACCTGTGTATATTTGAAGAAGTGTTGAAC GAGTGGATCCAGAGGGGACTGATAAGACGGTGCTGAAGTACAAGTGCCACACGATGAAGAAGCTGAGCAT GACAAGGACGCTTCTGACAGAAAAGAAGGAAGGAGAAGAAAACATAGATGGCGTAGAGTGGCTGCAAATG AAGGATAATGATTTCTCTATAGACCCAACATGATTTGCAGCTTTCTCCATGAACATGAAGATGAGGCTG CAGGGCCGGCCTCTGATAAACCTGTGGAGCTGGAAGACCAAGAGAGTCAGATGAAAGGCAGTGCGGAGCA GACACAGGAGCAGGAGAAGGTGGAGCACTCGGAAGTGGAGGACCCTGCTCCTTCAGGGGAGACCCCTTCT GAGATGGAAAGTTCTGTTTTTCATGGGAAGTCAAGATGGAAATGTTTCTCAAACAACCCAGTCGTGAAC TAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
Restriction Sites:	Sfil-Sfil
ACCN:	NM_026113
Insert Size:	684 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC048446](#), [AAH48446](#)

RefSeq Size: 729 bp

RefSeq ORF: 684 bp

Locus ID: 67371

UniProt ID: [Q9D8P7](#)

Cytogenetics: 10 B1

Gene Summary: Involved in RNA polymerase III-mediated transcription. Integral, tightly associated component of the DNA-binding TFIIIC2 subcomplex that directly binds tRNA and virus-associated RNA promoters (By similarity).[UniProtKB/Swiss-Prot Function]