

Product datasheet for MC205059

Cypt3 (NM_173367) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Cypt3
Synonyms:	1700006H02Rik; Ckt1r1
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC048527
 AGTCCCTGGCGAGCTCCGCTATGGCTCAGGTGGCCAAGAAAGTCCACGGGTCTAGAGCTGCGGCAGCAGC
 AGCAGCCGCGCTGCCAAAGCCAAAAGATCGAAGCCTAAGACCGCCGCAAAAGATTCAAGCTTAATAAG
 AAGCGAAACCCTAGATCTAAGCTTCCCAAGAGATCTCATCGTTCCTTAATTCGTTCTCTTTTCAGTCTC
 GTTCTGTGGTTGCTGTCGTTGTTGCTGTGGTTTTGCTGTTATCGTTCTCGTCCTTCTCGTAGCAGAAT
 TACATTTAAGATTACGAAGAAGGGGAACAATCACTTAGAAGGAGAATTAGGAGACGAATCAAGAGTGCT
 ACTGAACCTTCGATTGATGCGGAGCCAACTGGAGCTGATCGAACCGGAACCGACTATGGCTTTGGAGCCAA
 GTGAGATTACAGTGGCATTGTTCTCTCATGAGAATGCCAATGTGTCTGAGCCGGAGGAAGTCCACCATG
 CATTGACAGTGACCAGGTTCCGAATGGAGACTTGGCCAGCTCCTAGACCCTGAACAGGTCCAAAGACCCA
 TGAACATCAGTGAACGGGTCTGGGAGCTGAGCTGTGGTCTGTCAAGGCCACTTTTCACCACAACTAA
 TTTGTATGATGGTGAATAATAAGATGCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

Restriction Sites: Sfil-Sfil

ACCN: NM_173367

Insert Size: 516 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC048527](#), [AAH48527](#)

RefSeq Size: 689 bp

RefSeq ORF: 516 bp

Locus ID: 69361

Cytogenetics: X F3