

Product datasheet for **MC205154**

St6galnac5 (NM_012028) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	St6galnac5 (NM_012028) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	St6galnac5
Synonyms:	AI851940; Siat7e
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

RefSeq: [BC055737](#), [AAH55737](#)

RefSeq Size: 2026 bp

RefSeq ORF: 1008 bp

Locus ID: 26938

UniProt ID: [Q9QYJ1](#)

Cytogenetics: 3 H3

Gene Summary: Involved in the biosynthesis of ganglioside GD1a from GM1b. It exhibits higher activity with glycolipids than with glycoproteins.[UniProtKB/Swiss-Prot Function]