

Product datasheet for **MC205182**

Crygb (NM_144761) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Crygb (NM_144761) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Crygb
Synonyms:	Cryg-3; DGcry-3; Nop
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC056455 ACTCTCAGCGAGATGGGAAAGATCACCTTCTTCGAGGACCGCAGCTTCCAGGGCCGCTGCTATGAGTGCA GCAGCGACTGCCCAACCTGCAGACCTACTTCAGCCGCTGCAATTCTGTCCGCGTGGACAGTGGCTGCTG GATGCTCTATGAGCGCCCAACTACCAGGGCCACCACTACTTCCCTGAGACGTGGAGAGTACCCTGACTAC CAGCAGTGGATGGGTTTCAGCGACTCCATTTCGTTCTGCTGCCTCATCCCCAACACTCTGGCACTTACA GAATGAGGATCTACGAAAAAGATGACTTCAGAGGACAAATGTCAGAGATCACAGACGATTGTCTCTCTCT TCAGGATCGCTTCCACTTCAGTGAATCCACTCCCTCAATGTGATGGAGGGCTGCTGGGTTCTCTATGAG ATGCCTAGCTACAGAGGACGGCAGTACCTGCTAAGGCCGGGGGAGTACAGGAGATATCTTGACTGGGGGG CTGCAAATGCCAAAGTTGGCTCTTTTAGCGAGTCATGGATTTTACTGAGGCATTTTGAGACTCTACTT TCCTCCTCTAGAGCCAAATAAAGTATGTAGCTTATGATTCTGGCAAAAAAAAAAAAAAAAAAAAAAAAAAAAA AAAAAAAA
Restriction Sites:	Ascl-NotI
ACCN:	NM_144761
Insert Size:	528 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC056455](#), [AAH56455](#)

RefSeq Size: 638 bp

RefSeq ORF: 528 bp

Locus ID: 12965

UniProt ID: [P04344](#)

Cytogenetics: 1 32.82 cM

Gene Summary: Crystallins are the dominant structural components of the vertebrate eye lens.
[UniProtKB/Swiss-Prot Function]