

Product datasheet for MC205237

Sar1b (NM_025535) Mouse Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Sar1b (NM_025535) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Sar1b |
| Synonyms: | 2310075M17Rik; 2900019I22Rik; CMRD; Sara1b; Sara2; Sarb |
| Mammalian Cell Selection: | Neomycin |
| Vector: | PCMV6-Kan/Neo (PCMV6KN) |
| E. coli Selection: | Kanamycin (25 ug/mL) |

Fully Sequenced ORF: >BC082550
GAGAAGAGCTGGGGCCCTGGCAGCCATAGCGCTTTGCACAGCCGGCTTGGAGCGTGTAGGACGGGCACGG
AGCAGGGTTGCGTTCTCTGAAGTCGCAGACCTGTAAGATCTGGATATGTCCTTCATATTTGACTGGATTT
ACAGTGGTTTTAGCAGTGTGCTGCAGTTTTCTAGGATTGTATAAAAAGTCTGGTAAACTGGTATTCCTTGG
ATTGGATAATGCCGGGAAAAACAACCTTTGCTACACATGCTGAAAGATGACAGGCTCGGGCAGCACGTCCTCA
ACGCTACATCCCCTTCAGAAGAGCTTACTATTGCTGGCATGACGTTTACAACCTTTGATCTGGGTGGGC
ACGTGCAAGCCCGAAGAGTATGGAAAACTACCTTCTGCTATCAATGGCATTGTATTTCTGTGGATTG
TGCAGACCATGAAAGGCTGTTAGAATCAAAGAAGAACTTGATTCATAATGACAGATGAAACCATTGCC
AACGTGCCTATTCTAATTCTTGGAAATAAGATCGACAGACCTGAAGCCATCAGTGAAGAGAGGTTACGAG
AGATGTTTTGGTTTATATGGTCAGACAACAGGAAAGGGTAGTGTGTCAGTGAAGGAGCTGAACGCCGGCC
TCTGGAAGTGTTTATGTGCAAGTGTGCTAAAAAGGCAGGGCTATGGAGAAGGCTTCCGCTGGATGGCGCAG
TACATCGATTAGCAGCGCCGGCCACACACTTCCAGCCCTCAGAGATCCATCCCCAACATGCATAACT
TGAATTCATAGACTTTTGTGGTTTCAGAATAGGTGTTTTTAGATTATTAATATATCTACCTGAGTGGG
AATTGCAAAGTATTCAAGTTTGAATATCACAATGTTAGCTTTCTAATCCACACAGATAGTTTTTACAG
TTTCTAACCCACATCAATCCATCACCGTTTATAAAGATCAGCTTCCCATCAGTACATTTGAAGCACTTT
TTAACGATGTGAAACTACAAATATAAACCATATTTAAAAGCTCATCATGTGAAACTTTTTATGACTTTTT
GGAAATAGTTTTCAATTTTAGATTGTCTGCCTATCTTGTGTACAGTTAAAATCCGCTTATCAGTGATTG
CATGATGCCTTACGGTTTCCAAGACTACTTTTTCTTATGCAAACGTCATGCAATAAAAACAACCTCCAACA
TTTGAAAAAAAAAAAAAAAAAAAA

| | |
|--------------------|-----------|
| Restriction Sites: | Ascl-NotI |
| ACCN: | NM_025535 |
| Insert Size: | 597 bp |



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>BC082550</u> , <u>AAH82550</u> |
| RefSeq Size: | 1213 bp |
| RefSeq ORF: | 597 bp |
| Locus ID: | 66397 |
| UniProt ID: | <u>Q9CQC9</u> |
| Cytogenetics: | 11 B1.3 |
| Gene Summary: | Involved in transport from the endoplasmic reticulum to the Golgi apparatus. Activated by the guanine nucleotide exchange factor PREB. Involved in the selection of the protein cargo and the assembly of the COPII coat complex.[UniProtKB/Swiss-Prot Function] |