

Product datasheet for MC205241

Ltv1 (NM_181470) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ltv1 (NM_181470) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Ltv1
Synonyms:	2610020N02Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC067062

```

GTGACGTGCTGCCCTAGCCCTGTGGACCCACGTGTGGAGCTGCCGGCAGGTACCGGGGAGCGGGCGTTCCG
GCTGCCGCGGTTCCGGAGAGGGTGTCCCTGCCTCATCATGCCTCACAGAAAGAAAAGCCCTTTATAGAG
AAGAAGAAAGCTGTGTCTTTTCACTGGTCCACCGGAGCCAGCGAGATCCTTTAGCCGCAGATGAAACTG
CACCCAGAGGGTTCTCTTGCCACAAAAGGTCAATGACGAGGAAAGACGGGTGAGCAGAGGAAAGTA
CGGCGTGTCTTCGACGACGATTATGATTACTTGACGACCTGAAGGAGCCCTCTGGGCCCGCAGAGCTT
ATCCCCAGCACCTCCTTTGCCACACCCCTCTGCGACAAGACAGAAGACCCCTGCGTATACTCAAGCACTG
GAATTAAGTTGCCTTCATCTGTGTTTGCATCCGAGTTTGAGGAAGATGTTGGATTGTTAAACAAAGCAGC
TCCAGTTTCAGGACCTCGACTGGATTTTGTCCAGACATCGTTGCAGCTCTTGATGATGATTTTGTACTT
GATGACCCAGAAAACCTCCTTGAAGATGACTTTATTCTTCAGGCCAACAGCCGACAGGAGGAGAGAGGA
TGGATACGGAGTCTGAGGAAGACGATGGCCATGAGTGGGAAGACATGGACGACGAGGAAGGGAGTGATGA
CAGGAGCTCTGCAGGCTTCTGTGATGAGGTTGACCTATCTGCCCCCGGCTCACCGCAGGAAGCCATG
AAAAGACACCTGTTTTGGGAAGAGGAAACAAAAGCCGTTTTACAGAGTACTCCATGACCTCCTCAGTCA
TGAGGAGAAAATGAACAGCTGACTCTACACGATGAGAGGTTTGAGAAGTTTATGAACAGTATGATGATGT
TGAAATTTGAGCTCTGGATAATGCAGAATTGGAAGGTACCATTCAAGTGGACAGCAATCGGCTGCAGGAG
GTTTTGAATGACTACTACAAGGAGAAGGCAGAGAATTGTGTAATAATTGAGTACTCTTGAGCCTTTTGGAG
ATCAAGATCTGCCAACAAATGAGCTTGATGAGTCTGAAAAGGAAGAGACCATTACTGTAGTTCTTGAGGA
AGCCAAAGAGAAGTGGGATTGTGAATCCATTTGCAGTACATACTCCAATTTATACAACCACCCACAGCTC
ATCAAATATGAACCAAGCCCAACAAATCCACCTGTCGTCTAAGACCGGAATACCTCTCAATGTCTTAC
CTAAGAAAAGGCCTCACAGCAAAGCAGGTTGAACGGATGCAGATGATCAATGGCAGTGACCTACCCAAAGT
GTCCACCCAGCCTCGTTCCAAAGATGAGACCAAAGAGGATAAAAGAGCCAGGAAGCAAGCCATAAAAAGAG
GAGCGGAAGGAACGAAGAGTTGAAAAGAAAGCTAACAAAGCTAGCATTTAACTGGAGAAAAGAAGACAAG
AAAAGGAGCTACTGAACTTGAAGAAAAACATTGAAGGTCTAAAGCTGTAGGTGTGGGACATTCATATAAG
TTATATTTCCATGTCAGCAGTCTTAGTTATGCTCTGTAATTAACCAGGATCTGAAGAGGGGTGAGCTTGC
CCCTGCCATTTTTATATGCAAAATAGCTGAGAAAAATAGACTATTTGTATTTTCATACAAGCGTTTCTCT
GTCAGCTAACTTGTAATACTATAATGTTGTAAGTTAATATTACACGCTTATATTTCCCTGAAATAAA
TGATTTTATGAAAGTCAAAAAAAAAAAAAAAAAAAAA

```



[View online »](#)

Restriction Sites:	Ascl-NotI
ACCN:	NM_181470
Insert Size:	1413 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC067062</u> , <u>AAH67062</u>
RefSeq Size:	1784 bp
RefSeq ORF:	1413 bp
Locus ID:	353258
UniProt ID:	<u>Q6NSQ7</u>
Cytogenetics:	10 A2