

Product datasheet for **MC205303**

Pmel (NM_021882) Mouse Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Pmel (NM_021882) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Pmel |
| Synonyms: | D10H12S53E; D12S53Eh; gp87; gp100; Pmel17; Si; Silv |
| Mammalian Cell Selection: | Neomycin |
| Vector: | PCMV6-Kan/Neo (PCMV6KN) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

RefSeq: BC082555, AAH82555
RefSeq Size: 2039 bp
RefSeq ORF: 1881 bp
Locus ID: 20431
Cytogenetics: 10 77.13 cM