

## Product datasheet for **MC205343**

### Emc10 (NM\_197991) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Emc10 (NM\_197991) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Emc10  
**Synonyms:** 2310044H10Rik; 5430410O10Rik; Inm02; Mirta22  
**Mammalian Cell Selection:** Neomycin  
**Vector:** PCMV6-Kan/Neo (PCMV6KN)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >BC020179

```
GCGGGCTGTATGGCTCTCGGTTTTCTCAACGCTCCCCTATGGTGGCCGCGGGTGCCGGGGTGACCCGGC
TGTTAGTGCTCTTGCTGATGGTAGCCGCGCTCCTAGCAGAGCCGAGGCAGCGGCTGCCGGGTCGGGGC
CTCCGCGCGTGGGACCGGGCCGATGGCCGTGAAGCTGAGGGCTGTGGCACCGTGCTTTTCTGCTGGAG
CATTTCATTTGAGCTCGGTGATGGAGCCAACTTCCAGAAGCGAGGCTCGTCTGCTGGAACACAGCAGGATG
GCACCCTGTGGCAACACAGCGACAGCTCAGTGAGGAGGAGCGTGGCCGACTCCGGGATGTGGTCTGTG
CAATGGCCTCTACAGGGTCCGGTCCCAGGGCGGCTGGGACACTTGATGGTTCAGAAGCTGGCGGCCAT
GTGTCTTCTTCCGTCAGCGTCTCCCTGGTGGAGTCCGACCTTTCCGACCAGCTGACCTTGCACGTGG
ATGTGGCTGGCAACGTGGTGGCCTGTCTGTGGTGGTGTACCCTGGGGGCTGCCGGGGTCCGAGGTGGA
AGATGAGGACCTGGAGCTGTTCAATACATCTGTGCAGCTGCGGCCTCCAGCACTGCTCCAGGCCCGAG
ACTGCCGCTTTCATTGAGCGCCTGGAAATGGAGCAGGCCAGAAGGCTAAGAACCACAGGAGCAGAAGT
CTTTCTTTGCCAAATACTGGATGTACATCATTCCAGTTGTGCTGTTCCCTCATGATGTCAGGAGCGCCGGA
CGCCGGGGGCCAGGGTGGCGGTGGGGGGGGGAGCAGCCGGTGGAGCAGCTGTGCCACCTAGAACCCCC
CTCCCCCAGAGCCAGCCCAAGAAGGAGTTCTGTCCCACTCTCCCTATTGCATGAATATGGAAGGC
TGTCCCTTCAGTGAGCCCTCTGGCCTTCTGTAAAGCCCTCTTTCTGTCCCTGAGCCTCTCTCATCTCT
GTTGACTGAGAGCTTGGGTGGACCTCCCTGTAGCCAGCTCACTGGCTTCCCTGAGCACTGCAACTGTGTC
CCACCATGTGGCACTGTGCTCCTCTGTCTGCTAAACACCCACCAGCCTGCCCTACCCACCCACCATAC
ACTTTGGAACTTGCCAAGCTCTCTCCAGCCTCTGTGCCCTTGGCCCTGCAGGCCCGTGGCCCTCACT
GTCACCTCCAGCCCTTTGCCAAGGATCTGTGGCCAGAGGCCTCTGCTTTAGTGGCTAGGTGAGCCTC
CAGCCCACTGTCCAGGTGGCATGCTGTCTTCTTTGCCCCCTCTCTGGTGGCCCAAGATACCATGGTGAC
CTATCACTATCCTTTCTGCCTTTGGATGTCCAGCCTGGATCTGTCAACAGGAGAGGATTGTGGCCTCC
ACGTTAGTCTGTGAATGCACACTTCGAGTGACTTGTGTGAGGTTTTGAGAGCCGTTTTGCACTAGCTG
CTCGACAGCTGTGGCATGGCCGTGCTCTTGCACATGCGCCGCTGTGGCATGGGAATTGCTGTGCAGCC
TCAGCTGTGTTGATGGCTGCTGATTAACGTCCCTAAACAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
```

**Restriction Sites:** RsrII-NotI



[View online »](#)

<b>ACCN:</b>	NM_197991
<b>Insert Size:</b>	777 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC020179</a> , <a href="#">AAH20179</a>
<b>RefSeq Size:</b>	1609 bp
<b>RefSeq ORF:</b>	777 bp
<b>Locus ID:</b>	69683
<b>UniProt ID:</b>	<a href="#">Q3TAS6</a>
<b>Cytogenetics:</b>	7 B3
<b>Gene Summary:</b>	Promotes angiogenesis and tissue repair in the heart after myocardial infarction. Stimulates cardiac endothelial cell migration and outgrowth via the activation of p38 MAPK, PAK and MAPK2 signaling pathways.[UniProtKB/Swiss-Prot Function]