

Product datasheet for **MC205389**

Larp1b (NM_001040399) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Larp1b (NM_001040399) Mouse Untagged Clone
Tag: Tag Free
Symbol: Larp1b
Synonyms: 1700108L22Rik; 4933421B21Rik; AA690246; Larp2
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC025528

```

CCACGCGTCCGCCACGCGTCCGCGGACGCGTGGGGGCTAGTGACTTCAGTGATATGGCAAACCTGGCCAA
CACCAGGTGAATTGGTTAACACGGGGTCTCAGAGTGTTATCAACCAAGGAAATAAAAAGCCACAAATTAG
AAAAGAAAAGAGAAGAGATTGAAAAAGAAGCAACAGTGAGAGCAAAGAGAACCAGGAAAGCAAACTA
GATGGTCTACTGAGAACATCAGCGAGGATGAGGCTCAGTCAAGCAGCCAGAGGAAGAGAGCTAATAAGC
ACAAATGGGTACCACTCCACTTAGATGATGTGAGACCAGACAGCCAGGAAAGACCTGGGTCCCGAACAG
CTCTAGATGTCAACCTGAAGCAAATAAACCAACGCACAGTAATAGGAAAAGTGATACACGAAGTTGGAGA
CGGGAGAGAGAAAAGAGAGATGATCAAGATGAAGTGCCAGTGTGAGAAGTGAGGGTGGTACCATCCGAG
GTTCCACTAGAGGTGGGGACGAGGCCGGGGCCGAGGAAGAGGCCGAGGCCGAGGAAACCTCGATTGAA
CTTTGATTATTCATATGGTTATCGAGAGCCTGGGAAAAGGACTGAGCAGCCATTTTACAGAGCTTAAT
ACCAGCATGATGTATTACTATGATGATGGCACAGGTGTGCGTGTGTATCCTGTTGAAGAAACACTGCTTA
AAGAGTACATTAAGCGCCAAATTGAGTATTACTTCAGCACAGAGAACTTGAACGGGACTTCTTTCTTCG
AAGAAAAGATGGACGAGCAGGGGTTCCCTGCCATTTCCCTGATTGCTGGTTTTACCGTGTGCAGGCTCTC
ACCACAAACCTTAATCTCATCTTGGAGGCCTTAAGGATAGCACAGAAGTAGAAATTGTGGATGAGAAAA
TGAGGAAAAGATAGAACCAGAAAAATGGCCAATTCCAGGCCCTCCTCCACGAAATGTGCCACAAACAGA
CTTCTCTCAGTTGATTGATTGTCGGGAGTTCATACCTGGCCAGGCCTTTGGCTCCCATACAGTCAGGGTG
ATGATCTACTGATGCTCTGGATGTCATATGATCGGAATATACAGTGGAGTAGGTGACTCCCTGAAATAC
CATGTAATATGAGAGCAAATCGTGCTGCAGAAGTGAGAGTGACCATACTGTCTTCTGTGCAAGAATGTC
TACGCCACCCCTTGTGGAGAAGTTGGGCCAGCCTGCAGTGGAGAAAACCACTTTACGATAATGATG
TGAGACAGACCAGAGGACTGCGGTTACTGTAGGAGCCTCCTTTCTATGTCATTAGCTTAGGCTGCACG
GGCTACTTCATTATGGTGGGAAAGTTAATTGTTATCATTGGGGAATACTGCAAATTTTATCTTTTGT
TCCAGGCATTTATACTTTGCCAGATCAATAAAGCAGAAAAAGAGAGGATTGAGTAGCATTTCACAAAT
GTTAGGTCCTTTGCTATCTGACTTGCTCTTAAATAGTATACTTTGCATAACTCACATTTGTCATTTTT
CACATGTACTATAAGATGTATTAATGTGAATACTTAAGGGTAAGAAATTAAGTTAAATAAATAAGGGC
ATTTTTCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_001040399
Insert Size:	1008 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC025528 , AAH25528
RefSeq Size:	1713 bp
RefSeq ORF:	1008 bp
Locus ID:	214048
Cytogenetics:	3 B