

Product datasheet for MC205431

Lrrc29 (NM_177449) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lrrc29 (NM_177449) Mouse Untagged Clone
Tag: Tag Free
Symbol: Lrrc29
Synonyms: A630024J02Rik; FBL9; FBXL9
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC034290

```
CTGCGATGGCGGACTCGCTGCCTCTGGAGATGCTCACATATATCCTGAGCTTCTGCCTCTGTCTCAGATCA
GAAAGAGGCCTCCCTTGTGAGTCGGACCTGGTACTGTGCAGCCCAGAATGCCCTTCGGGAGGTAAGACCT
CACTTGGTCCATTCTTACTTATTCTCTATTCTCCTTCTCATGGCCCTGCTTATCTATCATCTGTCTATC
CATCCTTATATACTTACACATTCATTCTACTAGTCTCGGGGTGCAGAGCTTGGAGTGCACAGGTGAATAG
AACATCATCTTCTAAGGAGCACATTGTCCAGTAAAGGAAAAAACAAGTTAGCGTGGTACGGTGTGTTC
GTGCCATATCAATGAAAGTGTGAAACCCAGACAAGGACACCAAGACGAGACCATGCAGTGGCATCAGGGC
AGGCCTGAGAGATAGGATGGGAACAAGTTACTCTTTCTAGAACGGGCATTAAGGTATGGAAATCCAC
ATAACTGGATTAAGTTTTTGTGTTGGATGTGGTAGGAGAAGAGTCAAGAATGTGGTCAAGATGGCAAGG
TGGCAGATGTGCCAGGCAGCAGAGAGAGCTTGGGCTTTCTCCAGTTACTGCTGTGGGGCCGTAAGAAA
TGGTAAGACAGGGCTGGAGAGATGGCTCAGCAGTTAAGAGCACTAACTACACTTCCAGAGTTCTGAGTT
CAATCCCAGCAACCACATGGTGGCTCACAAACATCTGT AATGGGAGCCGATGCCCTCTGTGTCTGAAGA
GAGGACAGTGTATTCACATACATAAAAATAAATAATTTGAAAAAAGAAAAGGTAAAGACAGGCAAGGACA
GTGTGAAACAGTCTTAGAACTCCATACTCCAGAAATGGGTTAGAGGGTAAAGATGGAGGTTGGGGAGGCAG
CCCAGCCAGGAGGGCAACACAAACATCTGTGGTATTTGAAAGCAATAGCTTGTGTTAATTTTTATTTAAT
TAGTTTTGATTTAGTTATAGTGAAGTTAAGGGTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTGCT
CATTAGCTATTTGTGTTCTTTATCCACTTAAGGTTAGGTTAGTAAATTTGCTTGTATGAACCTTTTA
TATATGAAAAATAAACTCTTAGCTGGGCAGTGTGGTGCATGCCTTTAGTCCCAGCAGTCAGGAAGCAG
AGGAAGGCAGATCTCTGAGTTCAAGGCCAGCCTGGTCTGCAGAGTAAGTTCTGGATAGCCAGGGCTGCA
CAGAAAAAAATCCTGTACAAAACAAACAAACAAACAAACAAACAAACAAACAAAGAAACAAAACAAAAGATT
ATTCTCTATATATTGGTTAAAATACATTCTATTTGCTCTTTCATTTTATTTGTAGATAGCAGTGTTTAC
ATATGAGGTACAATATTAAGCTTCTGTGTTCTGTTTCATCATTGGTTACTTAAAGATCCCTTCCCCCAGC
GATGTCATCAATGCATGCTCATATTTTTATTCTTTTAAAGAAATATAGCTGATTTTTAGAGTAGTGATTGG
TTCAATAAAAAATAAAATTTCTGTTTTATGGTTTTAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
```

Restriction Sites: RsrII-NotI



[View online »](#)

ACCN:	NM_177449
Insert Size:	291 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC034290</u> , <u>AAH34290</u>
RefSeq Size:	1608 bp
RefSeq ORF:	291 bp
Locus ID:	234684
Cytogenetics:	8 D3