

Product datasheet for **MC205438**

Grwd1 (NM_153419) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Grwd1 (NM_153419) Mouse Untagged Clone
Tag: Tag Free
Symbol: Grwd1
Synonyms: A301; AI504353; Grwd; Wdr28
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >BC028896

```

GGCGCGCAAGGGCCGGAGTCGCACGTGTGAAACCGGAGAACCATGGAGGCTGAAACCTGCGACCCAGGA
ACCGAAGGCCCATCCCAGGTCTACCTGCCCGCAGAGGGCCGCACTGAGCGAGGGAGAGGAGCTGGTAA
TGGATGAAGAGGCTACGTGCTGTACCATCGTGCGAAACTGACTGATGATGCTTCGGATGCACAATCTC
CATGGCACCAGGCCGTACCCCTCAGAGGCGAGTATGACGATGAGGAAGATGAAGATGAAGAGGATGAGG
AGGAGCAGAAACCTCAGCTGGAACCTGGCCATGGTTCCTACTACGGTGGCATCAACCGAGTTCGGGTATC
TTGGCTGGGGGAGGAACCTGTGGCTGGAGTATGGTCAGAGAAGGGCCAAGTGGAGGTGTTCCGCACTGCGG
CGGCTCCTGCAGGTGGTGGATGACCCCCAGGCCCTGGCCATCTTCTCCGAGATGAGCAGGCCCAATAA
AACCCATCTTCTCCTTTGCTGGCCACATGGGAGAGGGCTTTGCTCTTGACTGGTCCCTCGGGTGCCGG
CCGCTTGCTGACTGGTACTGTGAGAAGAAGTCCACCTCTGGACACCCAGGAGGTGGTTCCTGGAAT
GTGGACCAGAGGCCCTTTGTTGGCCACACACGGTCTGTGGAGGATCTGCAGTGGTCAACCACTGAGGACA
CGGTGTTTGCCTCCTGCTGCTGATGCCTCTATTTCGCATCTGGGACATCCGAGCAGCCCCAGGCAAGGC
CTGCATGCTCACACAGCCACTGCCATGATGGAGATGTTAATGTCATCAGCTGGAGCCGAGGGAGCCC
TTCTGCTCAGCGGTGGGGATGATGGTGCCTCAAAGTCTGGGACCTGCGACAGTTAAGTCTGGCTCCC
CTGTGGCCACCTTCAAGCAACATATGGCCCCGTGACCTCTGTTGAATGGCACCCCTCAGGACAGTGGGGT
CTTTGCAGCCTCAGGTGCTGACAACAGATCACACAATGGGACCTGGCAGTTGAACGGGACCCAGAGTGC
GGGAGACAGAGACAGACCCTGGGCTGGCGCACTCCCCAGCAGTTGCTGTTTCGTGCACCAAGGTGAGA
CTGACCTGAAGGAGCTACACTGGCACCCACAGTCCCCGGGGTCTGATCAGCACCGCCCTGTCAGGCTT
CACCGTCTCCGCACCATCAGTGTGGAGGTGGCCAACACTTGGCCTCAAGACTTGGCTCCACATGGGAA
CTGGATCGGTTAACTTGGCAAGGATCATTAGGCCTGTAGACTGCACCTCAGGAAGGAAGACTTGTTC AAC
TGGCCTACGTGGAACAGAAGGAACGACTCTACTGACCTGGTATCTCTAAGGATTTGGATTGACCCTTGG
CTCTGCCTTACACCAAAGACTAGGTGCAGGCAAAGATCCTTGCCTGTCCCAGATACGGCCTCAGGATGAG
ATTTTTGTCAGTGTGACTCCTGGCTTGACCCCAAGGACCTCTGACCAGGGACCTAATTCTGGCTGTTGA
TCAACTTCAGGAACCTTAAGGGTTTTGCCATAGCTGGGCCAAGTCTGGGTCTGTGCTCACTCCAGTGGAA
TTAATTCTCTCTGCTTCTGGCCAGTTAGATTTGGTTGAGATCCAAGTAGGGAGCCATCTGGGGTGCC
CTCCTGTGCTGACCACTACTGTACTGGAGGAGGTCTGGGGCTTCTGCTGTAATAAAGGACCACAATAAG
GCCCCAGACACCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAA

```



[View online »](#)

Restriction Sites:	RsrII-NotI
ACCN:	NM_153419
Insert Size:	1029 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC028896 , AAH28896
RefSeq Size:	1831 bp
RefSeq ORF:	1029 bp
Locus ID:	101612
UniProt ID:	Q810D6
Cytogenetics:	7 B3
Gene Summary:	Histone binding-protein that regulates chromatin dynamics and minichromosome maintenance (MCM) loading at replication origins, possibly by promoting chromatin openness.[UniProtKB/Swiss-Prot Function]