

Product datasheet for **MC205446**

Ccl21b (NM_011335) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ccl21b (NM_011335) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ccl21b
Synonyms: 6CKBAC1; 6Ckine; A1159701; CKb9; Scya21; Scya21a; SLC; TCA4
Mammalian Cell Selection: Neomycin
Vector: PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >BC025974
 CCCACGCGTCCGCAACTTGC GGCTGTCCATCTCACCTACAGCTCTGGTCTCATCCTCAACTCAACCACAA
 TCATGGCTCAGATGACTCTGAGCCTCCTTAGCCTGGTCTCTGCATCCCCTGGACCAAGG
 CAGTGATGGAGGGGTCAGGACTGCTGCCTTAAGTACAGCCAGAAGAAAATCCCTACAGTATTGTCCGA
 GGCTATAGGAAGCAAGAACCAAGTTTAGGCTGTCCATCCCGGCAATCCTGTTCTCACCCCGGAAGCACT
 CTAAGCCTGAGCTATGTGCAAACCTGAGGAAGGCTGGGTGAGAACCTGATGCGCCGCTGGACGACC
 TCCAGCCCCAGGAAACAAGCCCGGCTGCAGGAAGAACCAGGGAACCTCTAAGTCTGGAAGAAAGGA
 AAGGGCTCCAAGGGCTGCAAGAGAACTGAACAGACACAGCCCTCAAGAGGATAGCCAGTAGCCCGCTG
 GAGCCCAGGAGATCCCCACGAACCTCAAGCTGGGTGGTTCACGGTCCAACCTCACAGGCAAAGGGGAGC
 TAGAAAACAGACTCAGGAGCCCAAGCAGCCACCTCATGCTGGCCTCCGTCACACCCCTTGCCCTGCTTC
 AACCATCACATCTGCACGGCCATCCCTTTCTTACCTGGCGGAGCTGCCTTCCCTGGGGTAGACCTAGAGA
 GTCAGAAGAAAGAGTGTCTCCCAGGGAATGAGGAAGGAGACAGCAGGACTGTCCCCTCTAGGAGGCTACT
 CAGGTCCCAGACCTGAACCTGCTCTCCATGGCGCCCTTCCCTTGTCTTGCACCTATGATTTATACCTA
 ACTGAATAAAAAGGTGATCCAGCCTCAA
 AA

Restriction Sites: RsrII-NotI
ACCN: NM_011335
Insert Size: 402 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC025974</u> , <u>AAH25974</u>
RefSeq Size:	969 bp
RefSeq ORF:	402 bp
Locus ID:	100042493
UniProt ID:	<u>P84443</u>
Cytogenetics:	4 22.12 cM