

Product datasheet for MC205447

Nfrkb (NM_172766) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nfrkb (NM_172766) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Nfrkb
Synonyms:	A530090G11Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC033595

```

CCACGCGTCCGCGCTTTTCGGAGCCACGGCCGACGAACCCCGGTGAGACCCGCGGACTCTGCTGAGG
GTGTATTTTATGTAGTCTACCTGTGGAGAAGAAGCAATGGATTCTTGGACCACATGCTGACAGATCCC
CTTGAACCTTGGTCTTGTGGAGATGGTACAGTACAGGAATCATGGAGGACTGCCTTCTGGGTGGTACCA
GGTTAGTCTGCCGAAGACCTTCTGGAGATCCCGAGATATTCTTTGATGTCGTCAGCCTGTCCACATG
GCAGGAAGTATTAAGTACTCTCAGCGTGAACCTCCAGCAGTTCCTGCCTCGTTCCTGCGGACAGT
GTGGAGCAGCAGAGGGAGCTCATCTTGCCTATTTCAGTGGGAGAACTCCGCTTCGGAACCCGCTGC
ACATCGCCGAGAAGCTCTCCGAGATGGACACTTAAACCTGAGGTGGTCAAGTATCGTCAGCTCTGCTT
CAAGTCACAGTACAAGCGGTATCTCAACTCCCAACAACAGTATTTCCACCGGCTGCTGAAACAGATTCTT
GCTTCCAGGAGTATCTTGGAGATGGCTCGTAGGAGTGGCCCTGCTCTTCCCTTCCACACAAGCACC
ATTCACCGTCCCAGCCGAGGAGCGGGAGTGGCGGACCCAGCAGCGTTACTTAAAAGTCTTGAGGGA
AGTGAAGGAGGAATGTGGGGACTGCTCTGCTCTGATGAGGAGGATCTCAGCTCATGGCTTCCAAGC
TCTCCAGCACGTTCTCCTAGTCTGCGGTGCCCTGAGGGTGGTCCACGCTTTCACCCACGGATATGA
AACTGCAGATAAAATAGAACTGGGGGACAGTGACCTGAAGTTAATGTTAAAGAAGCACCATGAGAAGCG
GAAACATCAACCAGATCACCCAGACCTTTTACAGGGGACCTGACCCCTCAGTGACATCATGACTCGAGTA
AATGCCGGCAGGAAGGGCTCCCTAGCAGCCTTGATGATTTGGCTGTTCTTAAAAAAGGTGAAGGAAA
AGGAGGAGAAGAAAAAGAAGAAAATAAAGTTGATTAAGTCGGAAGCAGAGGATCTGGCTGAGCCCTAAG
CAGTACTGAAGGGGTCCCAACTCTCTCACAGGCCCTCTCCCTGGCAATATCGTCTATCAAGGAGGAA
CCCCTGGAAGACATCAAGCCTTGCCTTGGAAATCAATGAAATATCTCCAGTTTTTCTCTTCTGTTAG
AGATCTTGCTGCTGGAAGTACAGGCTAGCCTTCCCTATGCTGGAGGATCGAGTTTTGGACTGGCAGTCTT
TCCAGCCAGCTCCCTCAACAGCTGGTTCTGCTGCCCAACTGGGCTGAGTTGGTGGTGCCTGCTCTG
CAGTATCTTGCTGGAGAAAGCCGAGCGTTCTTCTAGTTTCTCTCCATTTGTTGAATCAAGAGAAAA
CCCAGCAGTGGAAATGCTTGGTCAATCTCAAGATAATGAAAAGGAATTAGCTGCTCTTCCACCTGTG
GCTGGAACCAAGACCAGGCTTCTGTAAGGAAAATGAAGACAGCTCAGATGCCATGACACCTGTCCCT
CGAGTAAGAACTGACTATGTGGTTCGGCTAGCACAGGAGAGGAGAAACGGGTTTTTCAAGAACAGGAGC
GTTACAGGTATAGCCAACCCATAAGGCATTACCTTTTCGATGCATGGCTTTGAGTCTGTGGTGGGGCC
AGTGAAGGGCGTGGTGGACAAGGAGACCTCCCTCAACAAGGCTCGTGAGCATTACTGCTGCGCTCTGAC

```



[View online »](#)

Insert Size:	3891 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC033595</u> , <u>AAH33595</u>
RefSeq Size:	5178 bp
RefSeq ORF:	3891 bp
Locus ID:	235134
UniProt ID:	<u>Q6PIJ4</u>
Cytogenetics:	9 A4
Gene Summary:	Binds to the DNA consensus sequence 5'-GGGGAATCTCC-3'. [UniProtKB/Swiss-Prot Function]