

## Product datasheet for **MC205687**

### Uqcrh (NM\_025641) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Uqcrh (NM_025641) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Uqcrh
Synonyms:	2210416J04Rik; 2310021J10Rik; 2610041P16Rik
Mammalian Cell Selection:	Neomycin
Vector:	PCMV6-Kan/Neo (PCMV6KN)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>BC011388 GTACTACTCTGGTTGCGCTTGTGGATCCCAAACCGGCTCCCGACATGGGACTAGAGGACGAACGAAAGA TGCTCACTGGATCTGGAGACCCCAAGAGGAAGAAGAGGAAGAACTAGTGGACCCCTAACACAGTGAG AGAGCACTGTGAACAGCTGGAGAAGTGTGTAAAGGCCCGGAACGACTAGAGTTGTGTGATAATCGCGTG TCTTCTCGGTACAGACAGAAGAGGATTGTACAGAGGAGCTCTTTGACTTCTTGCATGCACGGGACCACT GTGTGGCCACAAGCTCTTTAAAACTTGAAGTAAATGTGCAGATTTCGTCTCTCAGCCTTGCTACTGG GAATCAGGAACAGTTCCTTGTGGTTCTGGACGTCGGTGTCTGATGGAGTGAGTTCTCGAGAACCTCACTG ACTTTGGCGGTAGCTTCTTCTGTGTGACTAGCAGTGACTTCATCTTAATAAACTGATCTGCAAACCTTA AAAAAAAAAAAAAAAAAA
Restriction Sites:	RsrII-NotI
ACCN:	NM_025641
Insert Size:	270 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [BC011388](#), [AAH11388](#)

**RefSeq Size:** 506 bp

**RefSeq ORF:** 270 bp

**Locus ID:** 66576

**UniProt ID:** [P99028](#)

**Cytogenetics:** 4 D1

**Gene Summary:** This is a component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is part of the mitochondrial respiratory chain. This protein may mediate formation of the complex between cytochromes c and c1 (By similarity).  
[UniProtKB/Swiss-Prot Function]